



Be Different
Be Better

Variety Choice in Torque
Offer optimized rotation
control in clinical
treatment

Mini ORG
Lingual
Anterior 0.018"

Digital 1.0 A·db
BONDING TARY
Orthodontist's
Good Assistant

Compensation angle
Especially suitable
for teeth that not
fully erupted

PRODUCT CATALOGUE



Alex Zhang

Chairman of RS Medico
Senior Engineer

ORTHODONTIC DEVOTION · 35 YEARS
ONLY FOR A SUPERIOR BRACKET
ORTHODONTIC SOUL · THE PIONEER OF
BRACKET INDUSTRIALIZATION IN CHINA
ORTHODONTIC GOAL · WORLD-CLASS
INNOVATIVE COMPANY IN DENTAL FIELD

RS MEDICO was established in 2015 by Alex Zhang who was the former general manager of 3M ORJ.

RS MEDICO is a professional orthodontic manufacturer focusing on high quality ortho products.


The first company in China to adopt CNC milling technology to produce ortho brackets with 24 hours automatic production.

With ISO13485, CE&FDA, our products have been sold worldwide to over 50 countries.

RS MEDICO has gained 9 innovation patents, and 29 utility model patents.

OEM and ODM are welcome!

Looking for the partners for such projects!



RS MEDICO

DEVELOPMENT HISTORY

RS MEDICO is in line with international high and new technology, and has invested in the construction of an intelligent, information-based, high-quality international advanced production line, and the intelligent management of all equipment has witnessed the development of China's medical products industry with reliable RS MEDICO products.

2015.12

Establish



*To be the first rank innovation
company for dental products the world*



2016

Gained FDA



2017

Gained CE



2018

Marketed SL
Bracket



2019

Marketed A·db
Bonding Tray



2020

Marketed Customized
Bracket



2021

Honored National
High-tech Enterprise



Content

Digital 1.0

08 Digital Customized Products

08-13

08

Bracket

14 Self-Ligating Bracket 4.0

14-51

20 Self-Ligating Bracket

26 Queen Series Bracket

30 Ricketts Bracket

32 Off-Set Bracket

34 MUT Bracket

38 Sunflower Bracket

42 Mini ORG Lingual

44 Sapphire

47 Ceramic

50 Ceramic SL

14

Buccal Tube

52 Self-Ligating Buccal Tube

52-64

54 Queen Series

56 Mini Series

58 Daisy Series Bondable Buccal Tubes

62 Daisy Series Weldable Buccal Tubes

52

Molar Band

65 Molar Band

65-71

68 Molar Bands Preweld with Buccal Tube

70 Lingual Attachment for Bands

65

Accessories

72	Face Mask	72-95
76	High Pull System	
78	Face Bow	
80	Spring	
85	Crimpable Hook	
89	Lingual Button&Lingual Cleat	
93	Lingual Retainer	

72

Elastic

96	Ligature Tie	96-101
98	Power Chain	
100	O-Ring	
101	Other Elastic	

96

Instruments

102	Glass Mirror	102-111
106	Bracket positioning gauge/Arch wire turret	
107	Typodont	
108	Force Gauge/Tweezers	
109	Bracket opener/Retractor	
111	Orthodontic BlackBackground Photo Image Contrast	

102

Pliers

112-124

112

Wires

126-132


126



BE DIFFERENT BE BETTER

Digital 1.0

A•db BONDING TRAY

Orthodontist's
Good Assistant






Advantages of A • db

A•db Procedure

Be Different Be Better



1
Clinic

- Scan tooth model or make oral scan
- Obtain digital tooth model in stl. file
- Send the file to RS MEDICO




2
RS Medico

- Design a customized tray with ortho software
- Set the brackets position
- Send the design to clinic for check and confirmation



3
Clinic
&
RS Medico

- 3D print the tray
- Insert jiggers with RS Self-ligating brackets



4
Clinic

- Bonding the brackets directly



Product Parameters (Digital ortho 1.0)

SMM2-340a	Queen Self-Ligating Bracket MBT 345 with hook,6/7buccal tube 0.022" (28pcs/kit),with bracket Bonding Tray,teeth model,open tool.
SMQ2-340a	Queen Self-Ligating Bracket Q standard 345 with hook,6/7buccal tube 0.022" (28pcs/kit),with bracket Bonding Tray,teeth model,open tool.
SMN2-340a	Queen Self-Ligating Bracket High Torque 345 with hook,6/7buccal tube 0.022" (28pcs/kit),with bracket Bonding Tray,teeth model,open tool.
SML2-340a	Queen Self-Ligating Bracket Low Torque 345 with hook,6/7buccal tube 0.022" (28pcs/kit),with bracket Bonding Tray,teeth model,open tool.


BE DIFFERENT BE BETTER

Digital Ortho 2.0

Goal Oriented Digital Orthodontics


1. Personalized orthodontic treatment
2. Optional customized brackets
3. Customized arch wires





Self-Ligating Bracket 4.0





**Tangential
Angle of Jaw**
Reduce occlusal
interference



Arch Shape
Round corner
more comfortable



**Special
Design Cap**
Easy open



Slim Design
Lower profile



**Horizontal
Pull Hook**
Matching
ortho screw

Product Code (with wires)

RS MEDICO · Product catalogue

MBT Torque	SQM2-340K	345 w/hook 0.022"	SQM1-340K	345 w/hook 0.018"
	SQM2-330K	3 w/hook 0.022"	SQM1-330K	3 w/hook 0.018"
	SQM2-300K	no hook 0.022"	SQM1-300K	no hook 0.018"
Q Standard Torque	SQQ2-340K	345 w/hook 0.022"	SQQ1-340K	345 w/hook 0.018"
	SQQ2-330K	3 w/hook 0.022"	SQQ1-330K	3 w/hook 0.018"
	SQQ2-300K	no hook 0.022"	SQQ1-300K	no hook 0.018"
High Torque	SQH2-340K	345 w/hook 0.022"	SQH1-340K	345 w/hook 0.018"
	SQH2-330K	3 w/hook 0.022"	SQH1-330K	3 w/hook 0.018"
	SQH2-300K	no hook 0.022"	SQH1-300K	no hook 0.018"
Low Torque	SQL2-340K	345 w/hook 0.022"	SQL1-340K	345 w/hook 0.018"
	SQL2-330K	3 w/hook 0.022"	SQL1-330K	3 w/hook 0.018"
	SQL2-300K	no hook 0.022"	SQL1-300K	no hook 0.018"


1. Sold together with 1st & 2nd Molar bondable single MBT buccal tube
2. Select 8 pieces NiTi wires, upper lower each, 0.014" / 0.016" / 0.016X0.022" / 0.018" / 0.022" / 0.025"
3. Buy 5 sets and get 1 tool free

give



Self-Ligating Bracket 4.0

Reduce occlusal interference



Product Parameters (MBT)

RS MEDICO · Product catalogue

MBT	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	-7°	-7°	-7°	+10°	+17°	+17°	+10°	-7°	-7°	-7°
	Tip	0	0	+8°	+8°	+4°	+4°	+8°	+8°	0	0
	In/out mm	0.75	0.75	0.6	0.8	0.8	0.8	0.8	0.6	0.75	0.75
	Width mm	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8
	0.022" no hook	SQM2-3015	SQM2-3014	SQM2-3013	SQM2-3012	SQM2-3011	SQM2-3021	SQM2-3022	SQM2-3023	SQM2-3024	SQM2-3025
	0.022" with hook	SQM2-3115	SQM2-3114	SQM2-3113	\	\	\	\	SQM2-3123	SQM2-3124	SQM2-3125
	0.018" no hook	SQM1-3015	SQM1-3014	SQM1-3013	SQM1-3012	SQM1-3011	SQM1-3021	SQM1-3022	SQM1-3023	SQM1-3024	SQM1-3025
	0.018" with hook	SQM1-3115	SQM1-3114	SQM1-3113	\	\	\	\	SQM1-3023	SQM1-3124	SQM1-3125
	Mandibular										
	Tooth	45	44	43	42	41	31	32	33	34	35
	Torque	-17°	-12°	-6°	-6°	-6°	-6°	-6°	-6°	-12°	-17°
	Tip	+2°	+2°	+3°	0	0	0	0	+3°	+2°	+2°
	In/out mm	0.85	0.75	0.6	0.75	0.75	0.75	0.75	0.6	0.75	0.85
	Width mm	2.8	2.8	2.8	2.7	2.7	2.7	2.7	2.8	2.8	2.8
	0.022" no hook	SQM2-3045	SQM2-3044	SQM2-3043	SQM2-3042	SQM2-3041	SQM2-3031	SQM2-3032	SQM2-3033	SQM2-3034	SQM2-3035
	0.022" with hook	SQM2-3145	SQM2-3144	SQM2-3143	\	\	\	\	SQM2-3133	SQM2-3134	SQM2-3135
	0.018" no hook	SQM1-3045	SQM1-3044	SQM1-3043	SQM1-3042	SQM1-3041	SQM1-3031	SQM1-3032	SQM1-3033	SQM1-3034	SQM1-3035
0.018" with hook	SQM1-3145	SQM1-3144	SQM1-3143	\	\	\	\	SQM1-3133	SQM1-3134	SQM1-3135	

Product Parameters (Q Standard) RS MEDICO · Product catalogue

Q Standard	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	-11°	-11°	+6°	+6°	+17°	+17°	+6°	+6°	-11°	-11°
	Tip	0	0	+8°	+8°	+4°	+4°	+8°	+8°	0	0
	In/out mm	0.75	0.75	0.6	0.8	0.8	0.8	0.8	0.6	0.75	0.75
	Width mm	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8
	0.022" no hook	SQQ2-3015	SQQ2-3014	SQQ2-3013	SQQ2-3012	SQQ2-3011	SQQ2-3021	SQQ2-3022	SQQ2-3023	SQQ2-3024	SQQ2-3025
	0.022" with hook	SQQ2-3115	SQQ2-3114	SQQ2-3113	\	\	\	\	SQQ2-3123	SQQ2-3124	SQQ2-3125
	0.018" no hook	SQQ1-3015	SQQ1-3014	SQQ1-3013	SQQ1-3012	SQQ1-3011	SQQ1-3021	SQQ1-3022	SQQ1-3023	SQQ1-3024	SQQ1-3025
	0.018" with hook	SQQ1-3115	SQQ1-3114	SQQ1-3113	\	\	\	\	SQQ1-3123	SQQ1-3124	SQQ1-3125
Mandibular											
Tooth	45	44	43	42	41	31	32	33	34	35	
Torque	-17°	-12°	+7°	-3°	-3°	-3°	-3°	+7°	-12°	-17°	
Tip	+2°	+2°	+3°	0	0	0	0	+3°	+2°	+2°	
In/out mm	0.85	0.75	0.6	0.75	0.75	0.75	0.75	0.6	0.75	0.85	
Width mm	2.8	2.8	2.8	2.7	2.7	2.7	2.7	2.8	2.8	2.8	
0.022" no hook	SQQ2-3045	SQQ2-3044	SQQ2-3043	SQQ2-3042	SQQ2-3041	SQQ2-3031	SQQ2-3032	SQQ2-3033	SQQ2-3034	SQQ2-3035	
0.022" with hook	SQQ2-3145	SQQ2-3144	SQQ2-3143	\	\	\	\	SQQ2-3133	SQQ2-3134	SQQ2-3135	
0.018" no hook	SQQ1-3045	SQQ1-3044	SQQ1-3043	SQQ1-3042	SQQ1-3041	SQQ1-3031	SQQ1-3031	SQQ1-3033	SQQ1-3034	SQQ1-3035	
0.018" with hook	SQQ1-3145	SQQ1-3144	SQQ1-3143	\	\	\	\	SQQ1-3133	SQQ1-3134	SQQ1-3135	

Product Parameters (Low Torque) RS MEDICO · Product catalogue

Low Torque	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	-11°	-11°	-9°	-5°	+2°	+2°	-5°	-9°	-11°	-11°
	Tip	0	0	+8°	+8°	+4°	+4°	+8°	+8°	0	0
	In/out mm	0.75	0.75	0.6	0.8	0.8	0.8	0.8	0.6	0.75	0.75
	Width mm	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8
	0.022" no hook	SQL2-3015	SQL2-3014	SQL2-3013	SQL2-3012	SQL2-3011	SQL2-3021	SQL2-3022	SQL2-3023	SQQ2-3024	SQQ2-3025
	0.022" with hook	SQL2-3115	SQL2-3114	SQL2-3113	\	\	\	\	SQL2-3123	SQQ2-3124	SQQ2-3125
	0.018" no hook	SQL1-3015	SQL1-3014	SQL1-3013	SQL1-3012	SQL1-3011	SQL1-3021	SQL1-3022	SQL1-3023	SQQ1-3024	SQQ1-3025
	0.018" with hook	SQL1-3115	SQL1-3114	SQL1-3113	\	\	\	\	SQL1-3123	SQQ1-3124	SQQ1-3125
	Mandibular										
	Tooth	45	44	43	42	41	31	32	33	34	35
	Torque	-17°	-12°	0	-11°	-11°	-11°	-11°	0	-12°	-17°
	Tip	+2°	+2°	+3°	0	0	0	0	+3°	+2°	+2°
	In/out mm	0.85	0.75	0.6	0.75	0.75	0.75	0.75	0.6	0.75	0.85
Width mm	2.8	2.8	2.8	2.7	2.7	2.7	2.7	2.8	2.8	2.8	
0.022" no hook	SQL2-3045	SQL2-3044	SQL2-3043	SQL2-3042	SQL2-3041	SQL2-3031	SQL2-3032	SQL2-3033	SQL2-3034	SQL2-3035	
0.022" with hook	SQL2-3145	SQL2-3144	SQL2-3143	\	\	\	\	SQL2-3133	SQL2-3134	SQL2-3135	
0.018" no hook	SQL1-3045	SQL1-3044	SQL1-3043	SQL1-3042	SQL1-3041	SQL1-3031	SQL1-3032	SQL1-3033	SQL1-3034	SQL1-3035	
0.018" with hook	SQL1-3145	SQL1-3144	SQL1-3143	\	\	\	\	SQL1-3133	SQL1-3134	SQL1-3135	

Product Parameters (High Torque)

RS MEDICO · Product catalogue

High Torque	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	-4°	-4°	+11°	+13°	+22°	+22°	+13°	+11°	-4°	-4°
	Tip	0	0	+8°	+8°	+4°	+4°	+8°	+8°	0	0
	In/out mm	0.75	0.75	0.6	0.8	0.8	0.8	0.8	0.6	0.75	0.75
	Width mm	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8
	0.022" no hook	SQH2-3015	SQH2-3014	SQH2-3013	SQH2-3012	SQH2-3011	SQH2-3021	SQH2-3022	SQH2-3023	SQH2-3024	SQH2-3025
	0.022" with hook	SQH2-3115	SQH2-3114	SQH2-3113	\	\	\	\	SQH2-3123	SQH2-3124	SQH2-3125
	0.018" no hook	SQH1-3015	SQH1-3014	SQH1-3013	SQH1-3012	SQH1-3011	SQH1-3021	SQH1-3022	SQH1-3023	SQH1-3024	SQH1-3025
	0.018" with hook	SQH1-3115	SQH1-3114	SQH1-3113	\	\	\	\	SQH1-3123	SQH1-3124	SQH1-3125
Mandibular											
Tooth	45	44	43	42	41	31	32	33	34	35	
Torque	-17°	-5°	+11°	+11°	+11°	+11°	+11°	+11°	-5°	-17°	
Tip	+2°	+2°	+3°	0	0	0	0	+3°	+2°	+2°	
In/out mm	0.85	0.75	0.6	0.75	0.75	0.75	0.75	0.6	0.75	0.85	
Width mm	2.8	2.8	2.8	2.7	2.7	2.7	2.7	2.8	2.8	2.8	
0.022" no hook	SQH2-3045	SQH2-3044	SQH2-3043	SQH2-3042	SQH2-3041	SQH2-3031	SQH2-3032	SQH2-3033	SQH2-3034	SQM2-3035	
0.022" with hook	SQH2-3145	SQH2-3144	SQH2-3143	\	\	\	\	SQH2-3133	SQH2-3134	SQM2-3135	
0.018" no hook	SQH1-3045	SQH1-3044	SQH1-3043	SQH1-3042	SQH1-3041	SQH1-3031	SQH1-3032	SQH1-3033	SQH1-3034	SQM1-3035	
0.018" with hook	SQH1-3145	SQH1-3144	SQH1-3143	\	\	\	\	SQH1-3133	SQH1-3134	SQM1-3135	



Self-ligating Bracket

CNC Milling Workmanship

Ensure all technical data in precision and smooth slot in low friction

Precise Carving Arrow

Easy open and close stable performance and long life use

Laser Arrow on Clip

Guide fast bonding to correct position

Four Wings Remained


Enough tie-wing area for better rotation control in final treatment

Various Choice in Torque


Offer optimized rotation control in clinical treatment



ORTHODONTIC BRACKET



- Keep the space for four wings ligation
- The groove is accurate and smooth
- High, medium and low multi torque selection



Product Parameters (MBT)

RS MEDICO · Product catalogue

MBT	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	-7°	-7°	-7°	+10°	+17°	+17°	+10°	-7°	-7°	-7°
	Tip	0	0	+8°	+8°	+4°	+4°	+8°	+8°	0	0
	In/out mm	0.75	0.75	0.59	0.89	0.8	0.8	0.89	0.59	0.75	0.75
	Width mm	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8
	0.022" no hook	SMM2-3015	SMM2-3014	SMM2-3013	SMM2-3012	SMM2-3011	SMM2-3021	SMM2-3022	SMM2-3023	SMM2-3024	SMM2-3025
	0.022" with hook	SMM2-3115	SMM2-3114	SMM2-3113	\	\	\	\	SMM2-3123	SMM2-3124	SMM2-3125
	0.018" no hook	SMM1-3015	SMM1-3014	SMM1-3013	SMM1-3012	SMM1-3011	SMM1-3021	SMM1-3022	SMM1-3023	SMM1-3024	SMM1-3025
	0.018" with hook	SMM1-3115	SMM1-3114	SMM1-3113	\	\	\	\	SMM1-3023	SMM1-3124	SMM1-3125
	Mandibular										
	Tooth	45	44	43	42	41	31	32	33	34	35
	Torque	-17°	-12°	-6°	-6°	-6°	-6°	-6°	-6°	-12°	-17°
	Tip	+2°	+2°	+3°	0	0	0	0	+3°	+2°	+2°
	In/out mm	0.74	0.65	0.6	0.95	0.95	0.95	0.95	0.6	0.65	0.74
	Width mm	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8
	0.022" no hook	SMM2-3045	SMM2-3044	SMM2-3043	SMM2-3042	SMM2-3041	SMM2-3031	SMM2-3032	SMM2-3033	SMM2-3034	SMM2-3035
	0.022" with hook	SMM2-3145	SMM2-3144	SMM2-3143	\	\	\	\	SMM2-3133	SMM2-3134	SMM2-3135
	0.018" no hook	SMM1-3045	SMM1-3044	SMM1-3043	SMM1-3042	SMM1-3041	SMM1-3031	SMM1-3032	SMM1-3033	SMM1-3034	SMM1-3035
	0.018" with hook	SMM1-3145	SMM1-3144	SMM1-3143	\	\	\	\	SMM1-3133	SMM1-3134	SMM1-3135

Product Parameters (Q Standard) RS MEDICO · Product catalogue

Q Standard	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	-11°	-11°	+6°	+6°	+17°	+17°	+6°	+6°	-11°	-11°
	Tip	0	0	+8°	+8°	+4°	+4°	+8°	+8°	0	0
	In/out mm	0.75	0.75	0.59	0.89	0.8	0.8	0.89	0.59	0.75	0.75
	Width mm	2.8	2.8	2.8	2.8	2.8	2.8	2.8	3.0	2.8	2.8
	0.022" no hook	SMQ2-3015	SMQ2-3014	SMQ2-3013	SMQ2-3012	SMQ2-3011	SMQ2-3021	SMQ2-3022	SMQ2-3023	SMQ2-3024	SMQ2-3025
	0.022" with hook	SMQ2-3115	SMQ2-3114	SMQ2-3113	\	\	\	\	SMQ2-3123	SMQ2-3124	SMQ2-3125
	0.018" no hook	SMQ1-3015	SMQ1-3014	SMQ1-3013	SMQ1-3012	SMQ1-3011	SMQ1-3021	SMQ1-3022	SMQ1-3023	SMQ1-3024	SMQ1-3025
	0.018" with hook	SMQ1-3115	SMQ1-3114	SMQ1-3113	\	\	\	\	SMQ1-3123	SMQ1-3124	SMQ1-3125
Mandibular											
Tooth	45	44	43	42	41	31	32	33	34	35	
Torque	-17°	-12°	+7°	-3°	-3°	-3°	-3°	+7°	-12°	-17°	
Tip	+2°	+2°	+3°	0	0	0	0	+3°	+2°	+2°	
In/out mm	0.74	0.65	0.6	0.95	0.95	0.95	0.95	0.6	0.65	0.74	
Width mm	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	
0.022" no hook	SMQ2-3045	SMQ2-3044	SMQ2-3043	SMQ2-3042	SMQ2-3041	SMQ2-3031	SMQ2-3032	SMQ2-3033	SMQ2-3034	SMQ2-3035	
0.022" with hook	SMQ2-3145	SMQ2-3144	SMQ2-3143	\	\	\	\	SMQ2-3133	SMQ2-3134	SMQ2-3135	
0.018" no hook	SMQ1-3045	SMQ1-3044	SMQ1-3043	SMQ1-3042	SMQ1-3041	SMQ1-3031	SMQ1-3031	SMQ1-3033	SMQ1-3034	SMQ1-3035	
0.018" with hook	SMQ1-3145	SMQ1-3144	SMQ1-3143	\	\	\	\	SMQ1-3133	SMQ1-3134	SMQ1-3135	


Product Parameters (High Torque) RS MEDICO · Product catalogue

High Torque	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	-4°	-4°	+11°	+13°	+22°	+22°	+13°	+11°	-4°	-4°
	Tip	0	0	+8°	+8°	+4°	+4°	+8°	+8°	0	0
	In/out mm	0.75	0.75	0.59	0.89	0.8	0.8	0.89	0.59	0.75	0.75
	Width mm	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8
	0.022" no hook	SMH2-3015	SMH2-3014	SMH2-3013	SMH2-3012	SMH2-3011	SMH2-3021	SMH2-3022	SMH2-3023	SMH2-3024	SMH2-3025
	0.022" with hook	SMH2-3115	SMH2-3114	SMH2-3113	\	\	\	\	SMH2-3123	SMH2-3124	SMH2-3125
	0.018" no hook	SMH1-3015	SMH1-3014	SMH1-3013	SMH1-3012	SMH1-3011	SMH1-3021	SMH1-3022	SMH1-3023	SMH1-3024	SMH1-3025
	0.018" with hook	SMH1-3115	SMH1-3114	SMH1-3113	\	\	\	\	SMH1-3123	SMH1-3124	SMH1-3125
	Mandibular										
	Tooth	45	44	43	42	41	31	32	33	34	35
	Torque	-17°	-5°	+11°	+11°	+11°	+11°	+11°	+11°	-5°	-17°
	Tip	+2°	+2°	+3°	0	0	0	0	+3°	+2°	+2°
	In/out mm	0.74	0.65	0.6	0.95	0.95	0.95	0.95	0.6	0.65	0.74
Width mm	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	
0.022" no hook	SMH2-3045	SMH2-3044	SMH2-3043	SMH2-3042	SMH2-3041	SMH2-3031	SMH2-3032	SMH2-3033	SMH2-3034	SMM2-3035	
0.022" with hook	SMH2-3145	SMH2-3144	SMH2-3143	\	\	\	\	SMH2-3133	SMH2-3134	SMM2-3135	
0.018" no hook	SMH1-3045	SMH1-3044	SMH1-3043	SMH1-3042	SMH1-3041	SMH1-3031	SMH1-3032	SMH1-3033	SMH1-3034	SMM1-3035	
0.018" with hook	SMH1-3145	SMH1-3144	SMH1-3143	\	\	\	\	SMH1-3133	SMH1-3134	SMM1-3135	

Product Parameters (Low Torque)

RS MEDICO · Product catalogue

Low Torque	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	-11°	-11°	-9°	-5°	+2°	+2°	-5°	-9°	-11°	-11°
	Tip	0	0	+8°	+8°	+4°	+4°	+8°	+8°	0	0
	In/out mm	0.75	0.75	0.59	0.89	0.8	0.8	0.89	0.59	0.75	0.75
	Width mm	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8
	0.022" no hook	SML2-3015	SML2-3014	SML2-3013	SML2-3012	SML2-3011	SML2-3021	SML2-3022	SML2-3023	SMQ2-3024	SMQ2-3025
	0.022" with hook	SML2-3115	SML2-3114	SML2-3113	\	\	\	\	SML2-3123	SMQ2-3124	SMQ2-3125
	0.018" no hook	SML1-3015	SML1-3014	SML1-3013	SML1-3012	SML1-3011	SML1-3021	SML1-3022	SML1-3023	SMQ1-3024	SMQ1-3025
	0.018" with hook	SML1-3115	SML1-3114	SML1-3113	\	\	\	\	SML1-3123	SMQ1-3124	SMQ1-3125
	Mandibular										
	Tooth	45	44	43	42	41	31	32	33	34	35
	Torque	-17°	-12°	0	-11°	-11°	-11°	-11°	0	-12°	-17°
	Tip	+2°	+2°	+3°	0	0	0	0	+3°	+2°	+2°
In/out mm	0.74	0.65	0.6	0.95	0.95	0.95	0.95	0.6	0.65	0.74	
Width mm	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	
0.022" no hook	SML2-3045	SML2-3044	SML2-3043	SML2-3042	SML2-3041	SML2-3031	SML2-3032	SML2-3033	SML2-3034	SML2-3035	
0.022" with hook	SML2-3145	SML2-3144	SML2-3143	\	\	\	\	SML2-3133	SML2-3134	SML2-3135	
0.018" no hook	SML1-3045	SML1-3044	SML1-3043	SML1-3042	SML1-3041	SML1-3031	SML1-3032	SML1-3033	SML1-3034	SML1-3035	
0.018" with hook	SML1-3145	SML1-3144	SML1-3143	\	\	\	\	SML1-3133	SML1-3134	SML1-3135	



Product Code (with wires & tool) RS MEDICO · Product catalogue

MBT Torque	SMM2-3407	345 w/hook 0.022"	SMM1-3407	345 w/hook 0.018"
	SMM2-3307	3 w/hook 0.022"	SMM1-3307	3 w/hook 0.018"
	SMM2-3007	no hook 0.022"	SMM1-3007	no hook 0.018"
Q Standard Torque	SMQ2-3407	345 w/hook 0.022"	SMQ1-3407	345 w/hook 0.018"
	SMQ2-3307	3 w/hook 0.022"	SMQ1-3307	3 w/hook 0.018"
	SMQ2-3007	no hook 0.022"	SMQ1-3007	no hook 0.018"
High Torque	SMH2-3407	345 w/hook 0.022"	SMH1-3407	345 w/hook 0.018"
	SMH2-3307	3 w/hook 0.022"	SMH1-3307	3 w/hook 0.018"
	SMH2-3007	no hook 0.022"	SMH1-3007	no hook 0.018"
Low Torque	SML2-3407	345 w/hook 0.022"	SML1-3407	345 w/hook 0.018"
	SML2-3307	3 w/hook 0.022"	SML1-3307	3 w/hook 0.018"
	SML2-3007	no hook 0.022"	SML1-3007	no hook 0.018"



1. Sold together with 1st & 2nd Molar bondable single MBT buccal tube
2. Select 8 pieces NiTi wires, upper lower each, 0.014 " / 0.016 " / 0.016X 0.022" / 0.018X 0.025"
3. Open tool




Product Code (no wire & tool) RS MEDICO · Product catalogue


MBT Torque	SMM2-3402	345 w/hook 0.022"	SMM1-3402	345 w/hook 0.018"
	SMM2-3302	3 w/hook 0.022"	SMM1-3302	3 w/hook 0.018"
	SMM2-3002	no hook 0.022"	SMM1-3002	no hook 0.018"
Q Standard Torque	SMQ2-3402	345 w/hook 0.022"	SMQ1-3402	345 w/hook 0.018"
	SMQ2-3302	3 w/hook 0.022"	SMQ1-3302	3 w/hook 0.018"
	SMQ2-3002	no hook 0.022"	SMQ1-3002	no hook 0.018"
High Torque	SMH2-3402	345 w/hook 0.022"	SMH1-3402	345 w/hook 0.018"
	SMH2-3302	3 w/hook 0.022"	SMH1-3302	3 w/hook 0.018"
	SMH2-3002	no hook 0.022"	SMH1-3002	no hook 0.018"
Low Torque	SML2-3402	345 w/hook 0.022"	SML1-3402	345 w/hook 0.018"
	SML2-3302	3 w/hook 0.022"	SML1-3302	3 w/hook 0.018"
	SML2-3002	no hook 0.022"	SML1-3002	no hook 0.018"



1. Sold together with 1st & 2nd Molar bondable single MBT buccal tube



Queen Series Bracket




1 High Comfort
Low profile make patient more comfortable

2 Rigorous Size
CNC slot ensure precision and less friction

3 Ample Space
Ample area under tie-wing for easy ligation

4 Tight Fit
Contoured base adapt the tooth perfectly



Product Parameters (MBT)

RS MEDICO · Product catalogue

MBT	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	-7°	-7°	-7°	+10°	+17°	+17°	+10°	-7°	-7°	-7°
	Tip	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
	In/out mm	0.74	0.74	0.57	0.9	0.67	0.67	0.9	0.57	0.74	0.74
	Width mm	2.9	2.9	3.0	2.7	3.4	3.4	2.7	3.0	2.9	2.9
	0.022" no hook	BQM2-2015	BQM2-2014	BQM2-2013	BQM2-2012	BQM2-2011	BQM2-2021	BQM2-2022	BQM2-2023	BQM2-2024	BQM2-2025
	0.022" with hook	BQM2-2115	BQM2-2114	BQM2-2113	/	/	/	/	BQM2-2123	BQM2-2124	BQM2-2125
	0.018" no hook	BQM1-2015	BQM1-2014	BQM1-2013	BQM1-2012	BQM1-2011	BQM1-2021	BQM1-2022	BQM1-2023	BQM1-2024	BQM1-2025
	0.018" with hook	BQM1-2115	BQM1-2114	BQM1-2113	/	/	/	/	BQM1-2123	BQM1-2124	BQM1-2125
Mandibular											
Tooth	45	44	43	42	41	31	32	33	34	35	
Torque	-17°	-12°	-6°	-6°	-6°	-6°	-6°	-6°	-6°	-12°	-17°
Tip	+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°	
In/out mm	0.9	0.76	0.51	1.14	1.14	1.14	1.14	0.51	0.76	0.9	
Width mm	2.9	2.9	3.0	2.4	2.4	2.4	2.4	3.0	2.9	2.9	
0.022" no hook	BQM2-2045	BQM2-2044	BQM2-2043	BQM2-2042	BQM2-2041	BQM2-2031	BQM2-2032	BQM2-2033	BQM2-2034	BQM2-2035	
0.022" with hook	BQM2-2145	BQM2-2144	BQM2-2143	/	/	/	/	BQM2-2133	BQM2-2134	BQM2-2135	
0.018" no hook	BQM1-2045	BQM1-2044	BQM1-2043	BQM1-2042	BQM1-2041	BQM1-2031	BQM1-2032	BQM1-2033	BQM1-2034	BQM1-2035	
0.018" with hook	BQM1-2145	BQM1-2144	BQM1-2143	/	/	/	/	BQM1-2133	BQM1-2134	BQM1-2135	



Product parameters (Roth) RS MEDICO · Product catalogue

Roth	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	-7°	-7°	-3°	+7°	+12°	+12°	+7°	-3°	-7°	-7°
	Tip	0°	0°	+10°	+9°	+5°	+5°	+9°	+10°	0°	0°
	In/out mm	0.74	0.74	0.55	0.89	0.67	0.67	0.89	0.55	0.74	0.74
	Width mm	2.9	2.9	3.0	2.7	3.4	3.4	2.7	3.0	2.9	2.9
	0.022" no hook	BQR2-2015	BQR2-2014	BQR2-2013	BQR2-2012	BQR2-2011	BQR2-2021	BQR2-2022	BQR2-2023	BQR2-2024	BQR2-2025
	0.022" with hook	BQR2-2115	BQR2-2114	BQR2-2113	/	/	/	/	BQR2-2123	BQR2-2124	BQR2-2125
	0.018" no hook	BQR1-2015	BQR1-2014	BQR1-2013	BQR1-2012	BQR1-2011	BQR1-2021	BQR1-2022	BQR1-2023	BQR1-2024	BQR1-2025
	0.018" with hook	BQR1-2115	BQR1-2114	BQR1-2113	/	/	/	/	BQR1-2123	BQR1-2124	BQR1-2125
	Mandibular										
	Tooth	45	44	43	42	41	31	32	33	34	35
	Torque	-22°	-17°	-7°	-1°	-1°	-1°	-1°	-7°	-17°	-22°
	Tip	+0°	0°	+6°	0°	0°	0°	0°	+6°	0°	+0°
	In/out mm	0.88	0.77	0.51	1.14	1.14	1.14	1.14	0.51	0.77	0.88
	Width mm	2.9	2.9	3.0	2.4	2.4	2.4	2.4	3.0	2.9	2.9
	0.022" no hook	BQR2-2045	BQR2-2044	BQR2-2043	BQR2-2042	BQR2-2041	BQR2-2031	BQR2-2032	BQR2-2033	BQR2-2034	BQR2-2035
	0.022" with hook	BQR2-2145	BQR2-2144	BQR2-2143	/	/	/	/	BQR2-2133	BQR2-2134	BQR2-2135
	0.018" no hook	BQR1-2045	BQR1-2044	BQR1-2043	BQR1-2042	BQR1-2041	BQR1-2031	BQR1-2032	BQR1-2033	BQR1-2034	BQR1-2035
	0.018" with hook	BQR1-2145	BQR1-2144	BQR1-2143	/	/	/	/	BQR1-2133	BQR1-2134	BQR1-2135

Product parameters (Edge-wise) RS MEDICO · Product catalogue

Edge-wise	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
	Tip	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
	In/out mm	0.6	0.6	0.6	0.65	0.65	0.65	0.65	0.6	0.6	0.6
	Width mm	3.0	3.0	3.0	2.7	3.4	3.4	2.7	3.0	3.0	3.0
	0.022" no hook	BQE2-2015	BQE2-2014	BQE2-2013	BQE2-2012	BQE2-2011	BQE2-2021	BQE2-2022	BQE2-2023	BQE2-2024	BQE2-2025
	0.022" with hook	BQE2-2115	BQE2-2114	BQE2-2113	/	/	/	/	BQE2-2123	BQE2-2124	BQE2-2125
	0.018" no hook	BQE1-2015	BQE1-2014	BQE1-2013	BQE1-2012	BQE1-2011	BQE1-2021	BQE1-2022	BQE1-2023	BQE1-2024	BQE1-2025
	0.018" with hook	BQE1-2115	BQE1-2114	BQE1-2113	/	/	/	/	BQE1-2123	BQE1-2124	BQE1-2125
Mandibular											
Tooth	45	44	43	42	41	31	32	33	34	35	
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Tip	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
In/out mm	0.6	0.6	0.6	0.65	0.65	0.65	0.65	0.6	0.6	0.6	0.6
Width mm	3.0	3.0	3.0	3.0	2.4	2.4	2.4	3.0	3.0	3.0	3.0
0.022" no hook	BQE2-2045	BQE2-2044	BQR2-2043	BQE2-2042	BQE2-2041	BQE2-2031	BQE2-2032	BQE2-2033	BQE2-2034	BQE2-2035	BQE2-2035
0.022" with hook	BQE2-2145	BQE2-2144	BQR2-2143	/	/	/	/	BQE2-2133	BQE2-2134	BQE2-2135	BQE2-2135
0.018" no hook	BQE1-2045	BQE1-2044	BQR1-2043	BQE1-2042	BQE1-2041	BQE1-2031	BQE1-2032	BQE1-2033	BQE1-2034	BQE1-2035	BQE1-2035
0.018" with hook	BQE1-2145	BQE1-2144	BQR1-2143	/	/	/	/	BQE1-2133	BQE1-2134	BQE1-2135	BQE1-2135

Product Code (MBT)

RS MEDICO · Product catalogue

MBT Bracket (Micro-etched base) Set:5-5	Prescription	0.022"	0.018"
	345 with hooks	BQM2-2401	BQM1-2401
	3 with hooks	BQM2-2301	BQM1-2301
	no hook	BQM2-2001	BQM1-2001

Product Code (Roth)

RS MEDICO · Product catalogue

Roth Bracket Set:5-5	Prescription	0.022"	0.018"
	345 with hooks	BQR2-2401	BQR1-2401
	3 with hooks	BQR2-2301	BQR1-2301
	no hook	BQR2-2001	BQR1-2001

Product Code (Edge-wise)

RS MEDICO · Product catalogue

Edge-wise Bracket Set:5-5	Prescription	0.022"	0.018"
	345 with hooks	BQE2-2401	BQE1-2401
	3 with hooks	BQE2-2301	BQE1-2301
	no hook	BQE2-2001	BQE1-2001

Package**Small and exquisite**

Preferred sheet material
Recyclable, environment-friendly
and pollution-free


**Convenient transportation**

Convenient packing
transportation Fast

**Stacking space**

No extra space waste

Ricketts Bracket



Sand blasting treatment of the bottom plate has high bonding strength




Accurate expression of groove data



Light force correction makes patients more comfortable

Product parameters

RS MEDICO · Product catalogue




Ricketts Bracket	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	0°	0°	+7°	+14°	+22°	+22°	+14°	+7°	0°	0°
	Tip	0°	0°	+8°	+10°	+4°	+4°	+10°	+8°	0°	0°
	In/out mm	1.1	0.6	1	0.9	0.9	0.9	0.9	1	0.6	1.1
	Width mm	2.5	2.5	2.7	2.5	3.2	3.2	2.5	2.7	2.5	2.5
	0.018" no hook	/	/	/	BQK1-2012	BQK1-2011	BQK1-2021	BQK1-2022	/	/	/
	0.018" with hook	BQK1-2115	BQK1-2114	BQK1-2113	/	/	/	/	BQK1-2123	BQK1-2124	BQK1-2125
	Mandibular										
	Tooth	45	44	43	42	41	31	32	33	34	35
	Torque	0°	0°	-7°	0°	0°	0°	0°	-7°	0°	0°
	Tip	0°	0°	+3°	0°	0°	0°	0°	+3°	0°	0°
	In/out mm	1.1	0.6	1	0.9	0.9	0.9	0.9	1	0.6	1.1
	Width mm	2.5	2.5	2.7	2.2	2.2	2.2	2.2	2.7	2.5	2.5
	0.018" no hook	/	/	/	BQK1-2042	BQK1-2041	BQK1-2031	BQK1-2032	/	/	/
	0.018" with hook	BQK1-2145	BQK1-2144	BQK1-2143	/	/	/	/	BQK1-2133	BQK1-2134	BQK1-2135

Product code

RS MEDICO · Product catalogue

Queen Series Ricketts Bracket (sand blasted) Set:5-5	
Prescription	0.022"
345 w/hooks	BQK1-2401

Off-Set Bracket



Compensation angle

Especially suitable for teeth that not fully erupted

Product Parameters (MBT)

RS MEDICO · Product catalogue

MBT								
Torque	Offset	In/Out mm	M/D Width mm	Tooth	W/hook 0.022"	W/hook 0.018"	no hook 0.022"	no hook 0.018"
-7°	0°	0.75	2.9	14	BQM2-2114a	BQM1-2114a	BQM2-2014a	BQM1-2014a
-7°	0°	0.75	2.9	24	BQM2-2124a	BQM1-2124a	BQM2-2024a	BQM1-2024a
-12°	+2°	0.75	2.9	34	BQM2-2134a	BQM1-2134a	BQM2-2034a	BQM1-2034a
-12°	+2°	0.75	2.9	44	BQM2-2144a	BQM1-2144a	BQM2-2044a	BQM1-2044a
-7°	0°	0.75	2.9	15	BQM2-2115a	BQM1-2115a	BQM2-2015a	BQM1-2015a
-7°	0°	0.75	2.9	25	BQM2-2125a	BQM1-2125a	BQM2-2025a	BQM1-2025a
-17°	+2°	0.75	2.9	35	BQM2-2135a	BQM1-2135a	BQM2-2035a	BQM1-2035a
-17°	+2°	0.75	2.9	45	BQM2-2145a	BQM1-2145a	BQM2-2045a	BQM1-2045a

Product Parameters (Roth)

RS MEDICO · Product catalogue

Roth								
Torque	Offset	In/Out mm	M/D Width mm	Tooth	W/hook 0.022"	W/hook 0.018"	no hook 0.022"	no hook 0.018"
-7°	0°	0.75	2.9	14	BQR2-2114a	BQR1-2114a	BQR2-2014a	BQR1-2014a
-7°	0°	0.75	2.9	24	BQR2-2124a	BQR1-2124a	BQR2-2024a	BQR1-2024a
-17°	0°	0.75	2.9	34	BQR2-2134a	BQR1-2134a	BQR2-2034a	BQR1-2034a
-17°	0°	0.75	2.9	44	BQR2-2144a	BQR1-2144a	BQR2-2044a	BQR1-2044a
-7°	0°	0.75	2.9	15	BQR2-2115a	BQR1-2115a	BQR2-2015a	BQR1-2015a
-7°	0°	0.75	2.9	25	BQR2-2125a	BQR1-2125a	BQR2-2025a	BQR1-2025a
-22°	0°	0.75	2.9	35	BQR2-2135a	BQR1-2135a	BQR2-2035a	BQR1-2035a
-22°	0°	0.75	2.9	45	BQR2-2145a	BQR1-2145a	BQR2-2045a	BQR1-2045a

Product Parameters (Edge-wise)

RS MEDICO · Product catalogue

Edge-wise								
Torque	Offset	In/Out mm	M/D Width mm	Tooth	W/hook 0.022"	W/hook 0.018"	no hook 0.022"	no hook 0.018"
0°	0°	0.60	2.9	14	BQE2-2114a	BQE1-2114a	BQE2-2014a	BQE1-2014a
0°	0°	0.60	2.9	24	BQE2-2124a	BQE1-2124a	BQE2-2024a	BQE1-2024a
0°	0°	0.60	2.9	34	BQE2-2134a	BQE1-2134a	BQE2-2034a	BQE1-2034a
0°	0°	0.60	2.9	44	BQE2-2144a	BQE1-2144a	BQE2-2044a	BQE1-2044a
0°	0°	0.60	2.9	15	BQE2-2115a	BQE1-2115a	BQE2-2015a	BQE1-2015a
0°	0°	0.60	2.9	25	BQE2-2125a	BQE1-2125a	BQE2-2025a	BQE1-2025a
0°	0°	0.60	2.9	35	BQE2-2135a	BQE1-2135a	BQE2-2035a	BQE1-2035a
0°	0°	0.60	2.9	45	BQE2-2145a	BQE1-2145a	BQE2-2045a	BQE1-2045a

MUT Brackets



Product Parameters (Roth)

RS MEDICO · Product catalogue

Roth

Maxillary							
	Torque	Ang.	IN/Out		Hook	L/R	0.22"
Central	+12°	5°	.034	0.86		UL	BQR4-8011
	+12°	5°	.034	0.86		UR	BQR4-8021
Lateral	+8°	9°	.044	1.12		UL	BQR4-8012
	+8°	9°	.044	1.12		UR	BQR4-8022
Cuspid	0°	11°	.032	0.81		UL	BQR4-8013
	0°	11°	.032	0.81		UR	BQR4-8023
	0°	11°	.032	0.81	D	UL	BQR4-8113
	0°	11°	.032	0.81	D	UR	BQR4-8123
Bicuspid	-7°	0°	.035	0.89		UNIV	BQR4-8014
	-7°	0°	.035	0.89	D	UL	BQR4-8114
	-7°	0°	.035	0.89	D	UR	BQR4-8124
Mandibular							
	Torque	Ang.	IN/Out		Hook	L/R	0.22"
Anterior	-1°	0°	.045	1.14		UNLV	BQR4-8031
Cuspid	-11°	5°	.030	0.76		LL	BQR4-8033
	-11°	5°	.030	0.76		LR	BQR4-8043
	-11°	5°	.030	0.76	D	LL	BQR4-8133
	-11°	5°	.030	0.76	D	LR	BQR4-8143
First Bicuspid	-17°	0°	.032	0.81		UNIV	BQR4-8034
	-17°	0°	.032	0.81	D	LL	BQR4-8134
	-17°	0°	.032	0.81	D	LR	BQR4-8144
Second Bicuspid	-22°	0°	.032	0.81		UNIV	BQR4-8035
	-22°	0°	.032	0.81	D	LL	BQR4-8135
	-22°	0°	.032	0.81	D	LR	BQR4-8145



Product Parameters (Edgewise) RS MEDICO · Product catalogue

Selection by Standard Edgewise Prescription	Maxillary							RS MEDICO · Product catalogue
		IN/Out					Adhesive Coated Appliance System	
		Torque	Ang.	In.	mm	Hook	L/R	0.22"
	Central	0°	0°	.030	0.76		UNIV	BQE4-8011
	Lateral	0°	0°	.030	0.76		UNIV	BQE4-8012
	Cuspid	0°	0°	.030	0.76		UNIV	BQE4-8013
		0°	0°	.030	0.76	D	UR	BQE4-8113
		0°	0°	.030	0.76	D	UL	BQE4-8123
	Bicuspid	0°	0°	.034	0.84		UNIV	BQE4-8014
		0°	0°	.034	0.84	D	UL/LR	BQE4-8114
0°		0°	.034	0.84	D	UR/LL	BQE4-8124	
Mandibular								
	IN/Out					Adhesive Coated Appliance System		
	Torque	Ang.	In.	mm	Hook	L/R	0.22"	
Anterior	0°	0°	.040	1.02		UNLV	BQE4-8031	
Cuspid	0°	0°	.030	0.76		LL	BQE4-8033	
	0°	0°	.030	0.76	D	LR	BQE4-8133	
	0°	0°	.030	0.76	D	LL	BQE4-8143	




Product Code (Roth) RS MEDICO · Product catalogue

Roth Bracket Set:5-5	Prescription	0.022"
	345 with hooks	BQR4-8401
	3 with hooks	BQR4-8301
	no hook	BQR4-8001

Product Code (Edge-wise) RS MEDICO · Product catalogue

Edge-wise Bracket Set:5-5	Prescription	0.022"
	345 with hooks	BQE4-8401
	3 with hooks	BQE4-8301
	no hook	BQE4-8001

Sunflower Bracket



Color-coded quadrant

Easy bracket identification



Laser Welding

Ensures more safety&healthy




Medium Size

Easy to control torque&comfortable



Torque-in-base design

Precise treatment



Product parameters (MBT)


RS MEDICO · Product catalogue

MBT	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	-7°	-7°	-7°	+10°	+17°	+17°	+10°	-7°	-7°	-7°
	Tip	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
	In/out mm	0.9	0.6	0.86	1.11	0.71	0.71	1.11	0.86	0.6	0.9
	Width mm	2.9	2.9	3.0	2.8	3.6	3.6	2.8	3.0	2.9	2.9
	0.022" no hook	BSM2-0015	BSM2-0014	BSM2-0013	BSM2-0012	BSM2-0011	BSM2-0021	BSM2-0022	BSM2-0023	BSM2-0024	BSM2-0025
	0.022" with hook	BSM2-0115	BSM2-0114	BSM2-0113	/	/	/	/	BSM2-0123	BSM2-0124	BSM2-0125
	0.018" no hook	BSM1-0015	BSM1-0014	BSM1-0013	BSM1-0012	BSM1-0011	BSM1-0021	BSM1-0022	BSM1-0023	BSM1-0024	BSM1-0025
	0.018" with hook	BSM1-0115	BSM1-0114	BSM1-0113	/	/	/	/	BSM1-0123	BSM1-0124	BSM1-0125
	Mandibular										
	Tooth	45	44	43	42	41	31	32	33	34	35
	Torque	-17°	-12°	-6°	-6°	-6°	-6°	-6°	-6°	-12°	-17°
	Tip	+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°
	In/out mm	0.63	0.6	0.85	1.1	1.1	1.1	1.1	0.85	0.6	0.63
	Width mm	2.9	2.9	3.0	2.4	2.4	2.4	2.4	3.0	2.9	2.9
	0.022" no hook	BSM2-0045	BSM2-0044	BSM2-0043	BSM2-0042	BSM2-0041	BSM2-0031	BSM2-0032	BSM2-0033	BSM2-0034	BSM2-0035
	0.022" with hook	BSM2-0145	BSM2-0144	BSM2-0143	/	/	/	/	BSM2-0133	BSM2-0134	BSM2-0135
	0.018" no hook	BSM1-0045	BSM1-0044	BSM1-0043	BSM1-0042	BSM1-0041	BSM1-0031	BSM1-0032	BSM1-0033	BSM1-0034	BSM1-0035
	0.018" with hook	BSM1-0145	BSM1-0144	BSM1-0143	/	/	/	/	BSM1-0133	BSM1-0134	BSM1-0135



Product parameters (Roth)

RS MEDICO · Product catalogue



Tooth position area
identification
Easy to identify

Roth	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	-7°	-7°	-3°	+7°	+12°	+12°	+7°	-3°	-7°	-7°
	Tip	0°	0°	+10°	+8°	+5°	+5°	+8°	+10°	0°	0°
	In/out mm	0.6	0.6	0.5	1.1	0.5	0.5	1.1	0.5	0.6	0.6
	Width mm	2.9	2.9	3.0	2.7	3.2	3.2	2.7	3.0	2.9	2.9
	0.022" no hook	BSR2-0015	BSR2-0014	BSR2-0013	BSR2-0012	BSR2-0011	BSR2-0021	BSR2-0022	BSR2-0023	BSR2-0024	BSR2-0025
	0.022" with hook	BSR2-0115	BSR2-0114	BSR2-0113	/	/	/	/	BSR2-0123	BSR2-0124	BSR2-0125
	0.018" no hook	BSR1-0015	BSR1-0014	BSR1-0013	BSR1-0012	BSR1-0011	BSR1-0021	BSR1-0022	BSR1-0023	BSR1-0024	BSR1-0025
	0.018" with hook	BSR1-0115	BSR1-0114	BSR1-0113	/	/	/	/	BSR1-0123	BSR1-0124	BSR1-0125
	Mandibular										
	Tooth	45	44	43	42	41	31	32	33	34	35
	Torque	-22°	-17°	-7°	-1°	-1°	-1°	-1°	-7°	-17°	-22°
	Tip	+0°	+0°	+6°	0°	0°	0°	0°	+6°	+0°	+0°
	In/out mm	0.5	0.5	0.5	0.9	0.9	0.9	0.9	0.5	0.5	0.5
	Width mm	2.85	2.85	3.0	2.4	2.4	2.4	2.4	3.0	2.85	2.85
	0.022" no hook	BSR2-0045	BSR2-0044	BSR2-0043	BSR2-0042	BSR2-0041	BSR2-0031	BSR2-0032	BSR2-0033	BSR2-0034	BSR2-0035
	0.022" with hook	BSR2-0145	BSR2-0144	BSR2-0143	/	/	/	/	BSR2-0133	BSR2-0134	BSR2-0135
	0.018" no hook	BSR1-0045	BSR1-0044	BSR1-0043	BSR1-0042	BSR1-0041	BSR1-0031	BSR1-0032	BSR1-0033	BSR1-0034	BSR1-0035
	0.018" with hook	BSR1-0145	BSR1-0144	BSR1-0143	/	/	/	/	BSR1-0133	BSR1-0134	BSR1-0135

Product parameters (Edge-wise)

RS MEDICO · Product catalogue



Medium size
comfortable

Edge-wise	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
	Tip	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
	In/out mm	0.6	0.6	0.6	0.7	0.7	0.7	0.7	0.6	0.6	0.6
	Width mm	3.0	3.0	3.0	2.7	3.4	3.4	2.7	3.0	3.0	3.0
	0.022" no hook	BSE2-0015	BSE2-0014	BSE2-0013	BSE2-0012	BSE2-0011	BSE2-0021	BSE2-0022	BSE2-0023	BSE2-0024	BSE2-0025
	0.022" with hook	BSE2-0115	BSE2-0114	BSE2-0113	/	/	/	/	BSE2-0123	BSE2-0124	BSE2-0125
	0.018" no hook	BSE1-0015	BSE1-0014	BSE1-0013	BSE1-0012	BSE1-0011	BSE1-0021	BSE1-0022	BSE1-0023	BSE1-0024	BSE1-0025
	0.018" with hook	BSE1-0115	BSE1-0114	BSE1-0113	/	/	/	/	BSE1-0123	BSE1-0124	BSE1-0125
	Mandibular										
	Tooth	45	44	43	42	41	31	32	33	34	35
	Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
	Tip	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
	In/out mm	0.6	0.6	0.6	0.65	0.65	0.65	0.65	0.6	0.6	0.6
	Width mm	3.0	3.0	3.0	2.4	2.4	2.4	2.4	3.0	3.0	3.0
	0.022" no hook	BSE2-0045	BSE2-0044	BSE2-0043	BSE2-0042	BSE2-0041	BSE2-0031	BSE2-0032	BSE2-0033	BSE2-0034	BSE2-0035
	0.022" with hook	BSE2-0145	BSE2-0144	BSE2-0143	/	/	/	/	BSE2-0133	BSE2-0134	BSE2-0135
	0.018" no hook	BSE1-0045	BSE1-0044	BSE1-0043	BSE1-0042	BSE1-0041	BSE1-0031	BSE1-0032	BSE1-0033	BSE1-0034	BSE1-0035
	0.018" with hook	BSE1-0145	BSE1-0144	BSE1-0143	/	/	/	/	BSE1-0133	BSE1-0134	BSE1-0135

Product code RS MEDICO · Product catalogue

MBT Bracket Set:5-5	Prescription	0.022"	0.018"
	345 with hooks	BSM2-0401	BSM1-0401
	3 with hooks	BSM2-0301	BSM1-0301
	no hook	BSM2-0001	BSM1-0001
Roth Bracket Set:5-5	Prescription	0.022"	0.018"
	345 with hooks	BSR2-0401	BSR1-0401
	3 with hooks	BSR2-0301	BSR1-0301
	no hook	BSR2-0001	BSR1-0001
Edge-wise Bracket Set:5-5	Prescription	0.022"	0.018"
	345 with hooks	BSE2-0401	BSE1-0401
	3 with hooks	BSE2-0301	BSE1-0301
	no hook	BSE2-0001	BSE1-0001



Mini ORG Lingual

Anterior 0.018" Posterior 0.022" design

Low profile

Low friction



Product code RS MEDICO · Product catalogue

Mini ORG Lingual Bracket

Low Profile and Low Friction

80 mesh base/Mini size/Anterior 0.018" slot, posterior 0.022" slot

tooth	11/12/13/21/22/23/31/32/33/41/42/43	44/45	16/26/17/27
slot size	0.018"	0.022"	0.022"

Product parameters RS MEDICO · Product catalogue

Mini ORG Lingual Bracket	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	-11°	-11°	-55°	-55°	-55°	-55°	-55°	-55°	-11°	-11°
	Tip	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
	In/out mm	0.7	0.7	1	1	1	1	1	1	0.7	0.7
	width mm	1.8	1.8	2.5	2.5	2.5	2.5	2.5	2.5	1.8	1.8
	0.022" no hook	/	/	TB71-2113	TB71-2111	TB71-2111	TB71-2111	TB71-2111	TB71-2113	/	/
	0.022" with hook	TB72-2114	TB72-2114	/	/	/	/	/	/	TB72-2114	TB72-2114
	Mandibular										
	Tooth	45	44	43	42	41	31	32	33	34	35
Torque	-11°	-11°	-40°	-40°	-40°	-40°	-40°	-40°	-40°	-11°	-11°
Tip	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
In/out mm	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7
width mm	1.8	1.8	2.2	2.2	2.2	2.2	2.2	2.2	2.2	1.8	1.8
0.022" no hook	/	/	TB71-2133	TB71-2131	TB71-2131	TB71-2131	TB71-2131	TB71-2133	/	/	
0.022" with hook	TB72-2114	TB72-2114	/	/	/	/	/	/	TB72-2114	TB72-2114	

Torque	Tip	In/out mm	Width mm	Tooth	0.022"
9°	0°	0.9	3.0	16	J672-8111
9°	0°	0.9	3.0	26	J672-8121
9°	0°	0.9	3.0	36	J672-8131
9°	0°	0.9	3.0	46	J672-8141
9°	0°	0.9	3.0	17	J772-5111
9°	0°	0.9	3.0	27	J772-5121
9°	0°	0.9	3.0	37	J772-5131
9°	0°	0.9	3.0	47	J772-5141

Set Code	Reference
TB73-24	5-5 20pcs/set
TB73-24J7	7-7 28pcs/set

Sapphire

Superior monocrystalline material ensures transparent and elegant

Smooth slot with lowest friction



Upgraded base design for easy bonding and debonding

Low contour&smooth surface provide patients with maximum comfort

Product parameters (MBT)

RS MEDICO · Product catalogue

MBT	Maxillary											
	Tooth	15	14	13	12	11	21	22	23	24	25	
	Torque	-7°	-7°	-7°	+10°	+17°	+17°	+10°	-7°	-7°	-7°	
	Tip	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°	
	In/out mm	1.0	1.0	1.1	1.2	1.0	1.0	1.2	1.1	1.0	1.0	
	Width mm	3.3	3.3	3.3	3.26	3.26	3.26	3.26	3.3	3.3	3.3	
	0.022" no hook	/	/	CHM2-0013	CHM2-0012	CHM2-0011	CHM2-0021	CHM2-0022	CHM2-0023	/	/	
	0.022" with hook	CHM2-0114	CHM2-0114	CHM2-0113	/	/	/	/	CHM2-0123	CHM2-0124	CHM2-0124	
	0.018" no hook	/	/	CHM1-0013	CHM1-0012	CHM1-0011	CHM1-0021	CHM1-0022	CHM1-0023	/	/	
	0.018" with hook	CHM1-0114	CHM1-0114	CHM1-0113	/	/	/	/	CHM1-0123	CHM1-0124	CHM1-0124	
	Mandibular											
	Tooth	45	44	43	42	41	31	32	33	34	35	
	Torque	-17°	-12°	-6°	-6°	-6°	-6°	-6°	-6°	-6°	-12°	-17°
	Tip	0°	0°	+3°	0°	0°	0°	0°	+3°	0°	0°	
	In/out mm	1.0	1.0	1.0	1.23	1.23	1.23	1.23	1.0	1.0	1.0	
	Width mm	3.3	3.3	3.3	2.88	2.88	2.88	2.88	3.3	3.3	3.3	
	0.022" no hook	/	/	CHM2-0043	CHM2-0031	CHM2-0031	CHM2-0031	CHM2-0031	CHM2-0033	/	/	
	0.022" with hook	CHM2-0145	CHM2-0144	CHM2-0143	/	/	/	/	CHM2-0133	CHM2-0134	CHM2-0135	
	0.018" no hook	/	/	CHM1-0043	CHM1-0031	CHM1-0031	CHM1-0031	CHM1-0031	CHM1-0033	/	/	
	0.018" with hook	CHM1-0145	CHM1-0144	CHM1-0143	/	/	/	/	CHM1-0133	CHM1-0134	CHM1-0135	

Product parameters (Roth)

RS MEDICO · Product catalogue

Roth	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	-7°	-7°	0°	+8°	+12°	+12°	+8°	0°	-7°	-7°
	Tip	0°	0°	+11°	+9°	+5°	+5°	+9°	+11°	0°	0°
	In/out mm	1.0	1.0	0.7	1.48	1.22	1.22	1.48	0.7	1.0	1.0
	Width mm	3.3	3.3	3.3	3.65	3.65	3.65	3.65	3.3	3.3	3.3
	0.022" no hook	/	/	CHR2-0013	CHR2-0012	CHR2-0011	CHR2-0021	CHR2-0022	CHR2-0023	/	/
	0.022" with hook	CHR2-0114	CHR2-0114	CHR2-0113	/	/	/	/	CHR2-0123	CHR2-0124	CHR2-0124
	0.018" no hook	/	/	CHR1-0013	CHR1-0012	CHR1-0011	CHR1-0021	CHR1-0022	CHR1-0023	/	/
	0.018" with hook	CHR1-0114	CHR1-0114	CHR1-0113	/	/	/	/	CHR1-0123	CHR1-0124	CHR1-0124
	Mandibular										
	Tooth	45	44	43	42	41	31	32	33	34	35
	Torque	-22°	-17°	-11°	-1°	-1°	-1°	-1°	-11°	-17°	-22°
	Tip	0°	0°	+5°	0°	0°	0°	0°	+5°	0°	0°
	In/out mm	1.0	1.0	0.82	1.23	1.23	1.23	1.23	0.82	1.0	1.0
	Width mm	3.05	3.05	3.3	2.88	2.88	2.88	2.88	3.3	3.05	3.05
	0.022" no hook	/	/	CHR2-0043	CHR2-0031	CHR2-0031	CHR2-0031	CHR2-0031	CHR2-0033	/	/
	0.022" with hook	CHR2-0145	CHR2-0144	CHR2-0143	/	/	/	/	CHR2-0133	CHR2-0134	CHR2-0135
	0.018" no hook	/	/	CHR1-0043	CHR1-0031	CHR1-0031	CHR1-0031	CHR1-0031	CHR1-0033	/	/
	0.018" with hook	CHR1-0145	CHR1-0144	CHR1-0143	/	/	/	/	CHR1-0133	CHR1-0134	CHR1-0135




- 1
Transparent cover designed
easy to recognize the brackets
- 2
Maxillary and mandibular
brackets packed separately,
easy to use in treatment

Product code RS MEDICO · Product catalogue

Sapphire Series MBT Brackets Set:5-5		
Prescription	0.022"	0.018"
345 with hooks	CHM2-0401	CHM1-0401
3 with hooks	CHM2-0301	CHM1-0301

Sapphire Series Roth Brackets Set:5-5		
Prescription	0.022"	0.018"
345 with hooks	CHR2-0401	CHR1-0401
3 with hooks	CHR2-0301	CHR1-0301

Ceramic



Product parameters (MBT) RS MEDICO · Product catalogue

MBT	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
	Tip	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
	In/out mm	0.7	0.7	0.7	1.1	0.7	0.7	1.1	0.7	0.7	0.7
	Width mm	3.0	3.3	3.6	3.3	3.7	3.7	3.3	3.6	3.3	3.0
	0.022" no hook	CMM2-0015	CMM2-0014	CMM2-0013	CMM2-0012	CMM2-0011	CMM2-0021	CMM2-0022	CMM2-0023	CMM2-0024	CMM2-0025
	0.022" with hook	CMM2-0115	CMM2-0114	CMM2-0113	/	/	/	/	CMM2-0123	CMM2-0124	CMM2-0125
	0.018" no hook	CMM1-0015	CMM1-0014	CMM1-0013	CMM1-0012	CMM1-0011	CMM1-0021	CMM1-0022	CMM1-0023	CMM2-0024	CMM2-0025
	0.018" with hook	CMM1-0115	CMM1-0114	CMM1-0113	/	/	/	/	CMM1-0123	CMM1-0124	CMM1-0124
	Mandibular										
	Tooth	45	44	43	42	41	31	32	33	34	35
	Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
	Tip	0°	0°	+3°	0°	0°	0°	0°	+3°	0°	0°
	In/out mm	0.7	0.7	0.7	1.1	1.1	1.1	1.1	0.7	0.7	0.7
	Width mm	3.3	3.3	3.4	3.1	3.1	3.1	3.1	3.4	3.3	3.3
	0.022" no hook	CMM2-0045	CMM2-0044	CMM2-0043	CMM2-0031	CMM2-0031	CMM2-0031	CMM2-0031	CMM2-0033	CMM2-0034	CMM2-0035
	0.022" with hook	CMM2-0145	CMM2-0144	CMM2-0143	/	/	/	/	CMM2-0133	CMM2-0134	CMM2-0135
	0.018" no hook	CMM1-0045	CMM1-0044	CMM1-0043	CMM1-0031	CMM1-0031	CMM1-0031	CMM1-0031	CMM1-0033	CMM1-0034	CMM1-0035
	0.018" with hook	CMM1-0145	CMM1-0144	CMM1-0143	/	/	/	/	CMM1-0133	CMM1-0134	CMM1-0135

Product parameters (Roth) RS MEDICO · Product catalogue

Roth	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	-7°	-7°	-2°	+7°	+12°	+12°	+7°	-2°	-7°	-7°
	Tip	0°	0°	+11°	+9°	+5°	+5°	+9°	+11°	0°	0°
	In/out mm	0.7	0.7	0.7	1.1	0.7	0.7	1.1	0.7	0.7	0.7
	Width mm	3.3	3.3	3.6	3.3	3.7	3.7	3.3	3.6	3.3	3.3
	0.022" no hook	CMR2-0015	CMR2-0014	CMR2-0013	CMR2-0012	CMR2-0011	CMR2-0021	CMR2-0022	CMR2-0023	CMR2-0024	CMR2-0025
	0.022" with hook	CMR2-0115	CMR2-0114	CMR2-0113	/	/	/	/	CMR2-0123	CMR2-0124	CMR2-0125
	0.018" no hook	CMR1-0015	CMR1-0014	CMR1-0013	CMR1-0012	CMR1-0011	CMR1-0021	CMR1-0022	CMR1-0023	CMR2-0024	CMR2-0025
	0.018" with hook	CMR1-0115	CMR1-0114	CMR1-0113	/	/	/	/	CMR1-0123	CMR1-0124	CMR1-0124
	Mandibular										
	Tooth	45	44	43	42	41	31	32	33	34	35
	Torque	-22°	-17°	-7°	-1°	-1°	-1°	-1°	-7°	-17°	-22°
	Tip	+0°	+0°	+6°	0°	0°	0°	0°	+6°	+0°	+0°
	In/out mm	0.7	0.7	0.7	1.1	1.1	1.1	1.1	0.7	0.7	0.7
	Width mm	3.3	3.3	3.4	3.1	3.1	3.1	3.1	3.4	3.3	3.3
	0.022" no hook	CMR2-0045	CMR2-0044	CMR2-0043	CMR2-0031	CMR2-0031	CMR2-0031	CMR2-0031	CMR2-0033	CMR2-0034	CMR2-0035
	0.022" with hook	CMR2-0145	CMR2-0144	CMR2-0143	/	/	/	/	CMR2-0133	CMR2-0134	CMR2-0135
	0.018" no hook	CMR1-0045	CMR1-0044	CMR1-0043	CMR1-0031	CMR1-0031	CMR1-0031	CMR1-0031	CMR1-0033	CMR1-0034	CMR1-0035
	0.018" with hook	CMR1-0145	CMR1-0144	CMR1-0143	/	/	/	/	CMR1-0133	CMR1-0134	CMR1-0135

Product code RS MEDICO · Product catalogue

Ceramic MBT Brackets Set:5-5		
Prescription	0.022"	0.018"
345 with hooks	CMM2-0401	CMM1-0401
3 with hooks	CMM2-0301	CMM1-0301
no hooks	CMM2-0001	CMM1-0001

Ceramic Roth Brackets Set:5-5		
Prescription	0.022"	0.018"
345 with hooks	CMR2-0401	CMR1-0401
3 with hooks	CMR2-0301	CMR1-0301



Product features

Ceramic SL

Smooth and Comfortable

Smooth appearance
lower profile, rounded shape

Transparent and Elegant

More transparent and invisible

New Base Design


Bonding strength is
greatly improved

New Material and New Technology

Imported high-purity polycrystalline
alumina and new CIM technology

Patented Locking System

Innovative passive sliding clip
easy operation



Product parameters (MBT)

RS MEDICO · Product catalogue

MBT	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
	Tip	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
	In/out mm	0.85	0.85	0.85	1.15	0.95	0.95	1.15	0.85	0.85	0.85
	width mm	3.2	3.2	3.2	3.2	3.5	3.5	3.2	3.2	3.2	3.2
	0.022" no hook	/	/	/	CSM2-0012	CSM2-0011	CSM2-0021	CSM2-0022	/	/	/
	0.022" with hook	CSM2-0115	CSM2-0114	CSM2-0113	/	/	/	/	CSM2-0123	CSM2-0124	CSM2-0125
	Mandibular										
	Tooth	45	44	43	42	41	31	32	33	34	35
	Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
	Tip	0°	0°	+3°	0°	0°	0°	0°	+3°	0°	0°
	In/out mm	0.8	0.8	0.85	1.15	1.15	1.15	1.15	0.85	0.8	0.8
	width mm	3.2	3.2	3.2	3.1	3.1	3.1	3.1	3.2	3.2	3.2
	0.022" no hook	/	/	/	CSM2-0031	CSM2-0031	CSM2-0031	CSM2-0031	CSM2-0033	/	/
	0.022" with hook	CSM2-0145	CSM2-0144	CSM2-0143	/	/	/	/	CSM2-0133	CSM2-0134	CSM2-0135

Product parameters (Roth)

RS MEDICO · Product catalogue

Roth	Maxillary										
	Tooth	15	14	13	12	11	21	22	23	24	25
	Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
	Tip	0°	0°	+11°	+9°	+5°	+5°	+9°	+11°	0°	0°
	In/out mm	0.85	0.85	0.85	1.15	0.95	0.95	1.15	0.85	0.85	0.85
	width mm	3.2	3.2	3.2	3.2	3.5	3.5	3.2	3.2	3.2	3.2
	0.022" no hook	/	/	CSR2-0013	CSR2-0012	CSR2-0011	CSR2-0021	CSR2-0022	CSR2-0023	/	/
	0.022" with hook	CSR2-0115	CSR2-0114	CSR2-0113	/	/	/	/	CSR2-0123	CSR2-0124	CSR2-0125
	Mandibular										
	Tooth	45	44	43	42	41	31	32	33	34	35
	Torque	-22°	-17°	-11°	-1°	-1°	-1°	-1°	-11°	-17°	-22°
	Tip	+0°	+0°	+5°	0°	0°	0°	0°	+5°	+0°	+0°
	In/out mm	0.8	0.8	0.85	1.15	1.15	1.15	1.15	0.85	0.8	0.8
	width mm	3.2	3.2	3.2	3.1	3.1	3.1	3.1	3.2	3.2	3.2
	0.022" no hook	/	/	/	CSR2-0031	CSR2-0031	CSR2-0031	CSR2-0031	/	/	/
	0.022" with hook	CSR2-0145	CSR2-0144	CSR2-0143	/	/	/	/	CSR2-0133	CSR2-0134	CSR2-0135

Ceramic Self-Ligating MBT Brackets Set: 5-5

Prescription	0.022"
345 with hooks	CSM2-0401
3 with hooks	CSM2-0301

Ceramic Self-ligating Roth Brackets Set: 5-5

Prescription	0.022"
345 with hooks	CSR2-0401
3 with hooks	CSR2-0301

Product features

Self Ligating Buccal Tube


Compact Overall Design



Central axis etching
line is easy to adhere
and accurate

Sand blasted mesh base


Cover stable and
reliable




Product parameters (MBT) RS MEDICO · Product catalogue

MBT	Torque	compensate	In/out mm	Width mm	Tooth	0.022"	0.022" (4pcs/kit)
	-14°	10°	0.65	3.8	16	JMM2-31163	JMM2-30163
	-14°	10°	0.65	3.8	26	JMM2-31263	
	-20°	0°	0.7	3.8	36	JMM2-31363	
	-20°	0°	0.7	3.8	46	JMM2-31463	

 Easy to use




 Small overall shape



Product features

Buccal Tube-Queen Series



MESIAL-DISTALVISIBLE ALIGNMENT LINE

For accurate replacing

DISTAL NOTCH


For easier holding

CONTOURED FUNNEL ENTRY

Wire easier passing

OCCCLUSAL INDENT

Fit better



Product parameters

RS MEDICO · Product catalogue



Queen Series 1st molar bondable single buccal tube									
	Torque	Offset	In/out mm	Width mm	Tooth	0.022"	0.018"	0.022" (4pcs/kit)	0.018" (4pcs/kit)
MBT	-14°	10°	0.65	3.8	16	TQM2-21160	TQM1-21160	TQM2-20160	TQM1-20160
	-14°	10°	0.65	3.8	26	TQM2-21260	TQM1-21260		
	-20°	0°	0.7	3.8	36	TQM2-21360	TQM1-21360		
	-20°	0°	0.7	3.8	46	TQM2-21460	TQM1-21460		
Roth	-10°	7°	0.65	3.8	16	TQR2-21160	TQR1-21160	TQR2-20160	TQR1-20160
	-10°	7°	0.65	3.8	26	TQR2-21260	TQR1-21260		
	-25°	6°	0.65	3.8	36	TQR2-21360	TQR1-21360		
	-25°	6°	0.65	3.8	46	TQR2-21460	TQR1-21460		
Edge-wise	0°	0°	0.6	3.8	16	TQE2-21160	TQE1-21160	TQE2-20160	TQE1-20160
	0°	0°	0.6	3.8	26	TQE2-21260	TQE1-21260		
	0°	0°	0.6	3.8	36	TQE2-21360	TQE1-21360		
	0°	0°	0.6	3.8	46	TQE2-21460	TQE1-21460		

BUCCAL TUBES

Product parameters

RS MEDICO · Product catalogue



Queen Series 2nd molar bondable single buccal tube									
	Torque	Offset	In/out mm	Width mm	Tooth	0.022"	0.018"	0.022" (4pcs/kit)	0.018" (4pcs/kit)
MBT	-14°	10°	0.65	3.2	17	TQM2-21170	TQM1-21170	TQM2-20170	TQM1-20170
	-14°	10°	0.65	3.2	27	TQM2-21270	TQM1-21270		
	-15°	0°	0.7	3.2	37	TQM2-21370	TQM1-21370		
	-15°	0°	0.7	3.2	47	TQM2-21470	TQM1-21470		
Roth	-10°	6°	0.65	3.2	17	TQR2-21170	TQR1-21170	TQR2-20170	TQR1-20170
	-10°	6°	0.65	3.2	27	TQR2-21270	TQR1-21270		
	-30°	6°	0.65	3.2	37	TQR2-21370	TQR1-21370		
	-30°	6°	0.65	3.2	47	TQR2-21470	TQR1-21470		
Edge-wise	0°	0°	0.6	3.2	17	TQE2-21170	TQE1-21170	TQE2-20170	TQE1-20170
	0°	0°	0.6	3.2	27	TQE2-21270	TQE1-21270		
	0°	0°	0.6	3.2	37	TQE2-21370	TQE1-21370		
	0°	0°	0.6	3.2	47	TQE2-21470	TQE1-21470		


Mini Buccal Tube

- **Ultra thin appearance**

More comfortable for patients

- **Distal clamping design**

Easier bonding of posterior teeth



- **Funnel shaped**

Guide the arch wire easily

- **Arc bottom plate design**

Especially suitable for teeth that without fully erupted



Product parameters (Mini Buccal tube)

RS MEDICO · Product catalogue

Mini Buccal tube									
	Torque	compensate	In/out mm	Width mm	Tooth	0.022"	0.018"	0.022" (4pcs/kit)	0.018" (4pcs/kit)
MBT	-14°	4°	0.7	2.5	17	TQM2-21170M	TQM1-21170M	TQM2-20170M	TQM1-20170M
	-14°	4°	0.7	2.5	27	TQM2-21270M	TQM1-21270M		
	-10°	0°	0.72	2.5	37	TQM2-21370M	TQM1-21370M		
	-10°	0°	0.72	2.5	47	TQM2-21470M	TQM1-21470M		
MBT-a	-14°	10°	0.7	2.5	17	TQM2-21170Ma	TQM1-21170Ma	TQM2-20170Ma	TQM1-20170Ma
	-14°	10°	0.7	2.5	27	TQM2-21270Ma	TQM1-21270Ma		
	-20°	0°	0.72	2.5	37	TQM2-21370Ma	TQM1-21370Ma		
	-20°	0°	0.72	2.5	47	TQM2-21470Ma	TQM1-21470Ma		
Roth	-14°	10°	0.7	2.5	17	TQR2-21170M	TQR1-21170M	TQR2-20170M	TQR1-20170M
	-14°	10°	0.7	2.5	27	TQR2-21270M	TQR1-21270M		
	-25°	4°	0.72	2.5	37	TQR2-21370M	TQR1-21370M		
	-25°	4°	0.72	2.5	47	TQR2-21470M	TQR1-21470M		

BUCCAL TUBES


Molar Bracket Product parameters

RS MEDICO · Product catalogue



Molar Bracket									
	Torque	Offset	In/out mm	Width mm	Tooth	0.022"	0.018"	0.022" (4pcs/kit)	0.018" (4pcs/kit)
MBT	-14°	10°	0.65	3.8	16	TQM2-31160	TQM1-31160	TQM2-30160	TQM1-30160
	-14°	10°	0.65	3.8	26	TQM2-31260	TQM1-31260		
	-20°	0°	0.65	3.8	36	TQM2-31360	TQM1-31360		
	-20°	0°	0.65	3.8	46	TQM2-31460	TQM1-31460		
Roth	-10°	7°	0.65	3.8	16	TQR2-31160	TQR1-31160	TQR2-30160	TQR1-30160
	-10°	7°	0.65	3.8	26	TQR2-31260	TQR1-31260		
	-25°	6°	0.65	3.8	36	TQR2-31360	TQR1-31360		
	-25°	6°	0.65	3.8	46	TQR2-31460	TQR1-31460		

Daisy Series Bondable Buccal Tubes



Product parameters (1st molar bondable single buccal tube) RS MEDICO · Product catalogue

1st molar bondable single buccal tube									
System	Torque	Offset	In/out mm	Width (M/D mm)	Tooth	Non-convertible		Convertible	
						0.022"	0.018"	0.022"	0.018"
MBT/Roth	-14°	10°	0.5	4.0	16	TDM2-21160	TDM1-21160	TDM2-11160	TDM1-11160
	-14°	10°	0.5	4.0	26	TDM2-21260	TDM1-21260	TDM2-11260	TDM1-11260
MBT	-20°	0°	0.5	4.0	36	TDM2-21360	TDM1-21360	TDM2-11360	TDM1-11360
	-20°	0°	0.5	4.0	46	TDM2-21460	TDM1-21460	TDM2-11460	TDM1-11460
Roth	-25°	4°	0.5	4.0	36	TDR2-21360	TDR1-21360	TDR2-11360	TDR1-11360
	-25°	4°	0.5	4.0	46	TDR2-21460	TDR1-21460	TDR2-11460	TDR1-11460
Edge-wise	0°	0°	0.5	4.0	16/36	TDE2-21160	TDE1-21160	TDE2-11160	TDE1-11160
	0°	0°	0.5	4.0	26/46	TDE2-21260	TDE1-21260	TDE2-11260	TDE1-11260




Convertible



Non-Convertible


Product parameters (1st molar bondable double rectangular buccal tube) RS MEDICO · Product catalogue

1st molar bondable double rectangular buccal tube									
System	Torque	Offset	In/out mm	Width (M/D mm)	Tooth	Non-convertible		Convertible	
						0.022"	0.018"	0.022"	0.018"
MBT/Roth	-14°	+10°	0.5	4	16	TDM2-22160	TDM1-22160	TDM2-12160	TDM1-12160
	-14°	+10°	0.5	4	26	TDM2-22260	TDM1-22260	TDM2-12360	TDM1-12360
MBT	-20°	0°	0.5	4	36	TDM2-22360	TDM1-22360	TDM2-12360	TDM1-12360
	-20°	0°	0.5	4	46	TDM2-22460	TDM1-22460	TDM2-12460	TDM1-12460
Roth	-25°	+4°	0.5	4	36	TDR2-22360	TDR1-22360	TDR2-12360	TDR1-12360
	-25°	+4°	0.5	4	46	TDR2-22460	TDR1-22460	TDR2-12460	TDR1-12460
Edge-wise	0°	0°	0.5	4	16/36	TDE2-22160	TDE1-22160	TDE2-12160	TDE1-12160
	0°	0°	0.5	4	26/46	TDE2-22260	TDE1-22260	TDE2-12260	TDE1-12260



Product parameters (1st molar bondable combination buccal tube "one round one rectangular")

1st molar bondable double buccal tube (one round one rectangular)									
System	Torque	Offset	In/out mm	Width (M/D mm)	Tooth	Non-convertible		Convertible	
						0.022"	0.018"	0.022"	0.018"
MBT/Roth	-14°	+10°	0.5	4	16	TDM2-23160	TDM1-23160	TDM2-13160	TDM1-13160
	-14°	+10°	0.5	4	26	TDM2-23260	TDM1-23260	TDM2-13360	TDM1-13360
MBT	-20°	0°	0.5	4	36	TDM2-23360	TDM1-23360	TDM2-13360	TDM1-13360
	-20°	0°	0.5	4	46	TDM2-23460	TDM1-23460	TDM2-13460	TDM1-13460
Roth	-25°	+4°	0.5	4	36	TDR2-23360	TDR1-23360	TDR2-13360	TDR1-13360
	-25°	+4°	0.5	4	46	TDR2-23460	TDR1-23460	TDR2-13460	TDR1-13460
Edge-wise	0°	0°	0.5	4	16/36	TDE2-23160	TDE1-23160	TDE2-13160	TDE1-13160
	0°	0°	0.5	4	26/46	TDE2-23260	TDE1-23260	TDE2-13260	TDE1-13260



BUCCAL TUBES

Product parameters (1st molar bondable triple buccal tube)

RS MEDICO · Product catalogue

1st molar bondable triple buccal tube									
System	Torque	Offset	In/out mm	Width (M/D mm)	Tooth	Non-Convertible		Convertible	
						0.022"	0.018"	0.022"	0.018"
MBT	-14°	10°	0.5	4.0	16	TDM2-24160	TDM1-24160	TDM2-14160	TDM1-14160
	-14°	10°	0.5	4.0	26	TDM2-24260	TDM1-24260	TDM2-14260	TDM1-14260
Roth	-14°	10°	0.5	4.0	16	TDR2-24160	TDR1-24160	TDR2-14160	TDR1-14160
	-14°	10°	0.5	4.0	26	TDR2-24260	TDR1-24260	TDR2-14260	TDR1-14260
Edge-wise	0°	0°	0.5	4.0	16	TDE2-24160	TDE1-24160	TDE2-14160	TDE1-14160
	0°	0°	0.5	4.0	26	TDE2-24260	TDE1-24260	TDE2-14260	TDE1-14260



Product parameters (Double Rectangular(Twin))

RS MEDICO · Product catalogue

Double Rectangular(Twin)							
System	Torque	Offset	In/out mm	Width (M/D mm)	Tooth	Non-Convertible	
						0.022"	0.018"
Roth	-14°	+10°	0.5	4.0	16	TDR2-22160e	TDR1-22160e
	-14°	+10°	0.5	4.0	26	TDR2-22260e	TDR1-22260e
Roth	-25°	+4°	0.5	4.0	36	TDR2-22360e	TDR1-22360e
	-25°	+4°	0.5	4.0	46	TDR2-22460e	TDR1-22460e
Edge-wise	0°	0°	0.5	4.0	16/36	TDE2-22160e	TDE1-22160e
	0°	0°	0.5	4.0	26/46	TDE2-22260e	TDE1-22260e




Buccal Tube (Bondable) in kit 4pcs/kit

RS MEDICO · Product catalogue

Buccal Tube (Bondable) in kit 4pcs/kit					
System	Kit Description	Convertible		Non-convertible	
		0.022"	0.018"	0.022"	0.018"
MBT	Upper Single Lower Single	TDM2-10160	TDM1-10160	TDM2-20160	TDM1-20160
	Upper Double Rectanglar Lower Single	TDM2-10260	TDM1-10260	TDM2-20260	TDM1-20260
	Upper One Round One Rectangular Lower Single	TDM2-10360	TDM1-10360	TDM2-20360	TDM1-20360
	Upper Triple Lower Single	TDM2-10460	TDM1-10460	TDM2-20460	TDM1-20460
	Upper Triple Lower Double Rectanglar	TDM2-10560	TDM1-10560	TDM2-20560	TDM1-20560
	Upper Triple Lower One Rectangular	TDM2-10660	TDM1-10660	TDM2-20660	TDM1-20660
	Upper Double Rectanglar Lower Double Rectanglar	TDM2-10760	TDM1-10760	TDM2-20760	TDM1-20760
Roth	Upper Single Lower Single	TDR2-10160	TDR1-10160	TDR2-20160	TDR1-20160
	Upper Double Rectanglar Lower Single	TDR2-10260	TDR1-10260	TDR2-20260	TDR1-20260
	Upper One Round One Rectangular Lower Single	TDR2-10360	TDR1-10360	TDR2-20360	TDR1-20360
	Upper Triple Lower Single	TDR2-10460	TDR1-10460	TDR2-20460	TDR1-20460
	Upper Triple Lower Double Rectanglar	TDR2-10560	TDR1-10560	TDR2-20560	TDR1-20560
	Upper Triple Lower One Rectangular	TDR2-10660	TDR1-10660	TDR2-20660	TDR1-20660
	Upper Double Rectanglar Lower Double Rectanglar	TDR2-10760	TDR1-10760	TDR2-20760	TDR1-20760
Edge-wise	Upper Single Lower Single	TDE2-10160	TDE1-10160	TDE2-20160	TDE1-20160
	Upper Double Rectanglar Lower Single	TDE2-10260	TDE1-10260	TDE2-20260	TDE1-20260
	Upper One Round One Rectangular Lower Single	TDE2-10360	TDE1-10360	TDE2-20360	TDE1-20360
	Upper Triple Lower Single	TDE2-10460	TDE1-10460	TDE2-20460	TDE1-20460
	Upper Triple Lower Double Rectanglar	TDE2-10560	TDE1-10560	TDE2-20560	TDE1-20560
	Upper Triple Lower One Rectangular	TDE2-10660	TDE1-10660	TDE2-20660	TDE1-20660
	Upper Double Rectanglar Lower Double Rectanglar	TDE2-10760	TDE1-10760	TDE2-20760	TDE1-20760




Daisy Series Weldable Buccal Tubes



Product parameters (1st Molar Double Rectangular weldable)

RS MEDICO · Product catalogue

1st Molar Double Rectangular Weldable									
Rx	Torque	Offset	In/out mm	Width (M/D mm)	R/L	Non-convertible		Convertible	
						0.022"	0.018"	0.022"	0.018"
Roth/MBT	-14°	+10°	0.5	4	UR	TDM2-62160	TDM1-62160	TDM2-52160	TDM1-52160
					UL	TDM2-62260	TDM1-62260	TDM2-52360	TDM1-52360
Roth	-25°	+4°	0.5	4	LL	TDR2-62360	TDR1-62360	TDR2-52360	TDR1-52360
					LR	TDR2-62460	TDR1-62460	TDR2-52460	TDR1-52460
MBT	-20°	0	0.5	4	LL	TDM2-62360	TDM1-62360	TDM2-52360	TDM1-52360
					LR	TDM2-62460	TDM1-62460	TDM2-52460	TDM1-52460
Edgewise	0	0	0.5	4	UR/LL	TDE2-62160	TDE1-62160	TDE2-52160	TDE1-52160
					UL/LR	TDE2-62260	TDE1-62260	TDE2-52260	TDE1-52260



Product parameters (1st molar single weldable)

RS MEDICO · Product catalogue


1st Molar Single Weldable									
Rx	Torque	Offset	In/out mm	Width (M/D mm)	R/L	Non-convertible		Convertible	
						0.022"	0.018"	0.022"	0.018"
Roth/MBT	-14°	+10°	0.5	4	UR	TDM2-61160	TDM1-61160	TDM2-51160	TDM1-51160
					UL	TDM2-61260	TDM1-61260	TDM2-51360	TDM1-51360
Roth	-25°	+4°	0.5	4	LL	TDR2-61360	TDR1-61360	TDR2-51360	TDR1-51360
					LR	TDR2-61460	TDR1-61460	TDR2-51460	TDR1-51460
MBT	-20°	0	0.5	4	LL	TDM2-61360	TDM1-61360	TDM2-51360	TDM1-51360
					LR	TDM2-61460	TDM1-61460	TDM2-51460	TDM1-51460
Edgewise	0	0	0.5	4	UR/LL	TDE2-61160	TDE1-61160	TDE2-51160	TDE1-51160
					UL/LR	TDE2-61260	TDE1-61260	TDE2-51260	TDE1-51260



Product parameters (2nd molar single weldable)


RS MEDICO · Product catalogue

2nd Molar Single Weldable								
Rx	Torque	Offset	In/out mm	Width (M/D mm)	R/L	Non-convertible		
						0.022"	0.018"	
Roth/MBT	-14°	+10°	0.5	4	UR	TDM2-61170	TDM1-61170	
					UL	TDM2-61270	TDM1-61270	
Roth	-25°	+4°	0.5	4	LL	TDR2-61370	TDR1-61370	
					LR	TDR2-61470	TDR1-61470	
MBT	-20°	0	0.5	4	LL	TDM2-61370	TDM1-61370	
					LR	TDM2-61470	TDM1-61470	
Edgewise	0	0	0.5	4	UR/LL	TDE2-61170	TDE1-61170	
					UL/LR	TDE2-61270	TDE1-61270	



Product parameters (1st molar Double Rectangular(Twin) Weldable) RS MEDICO · Product catalogue


1st Molar Double Rectangular(Twin) Weldable							
Rx	Torque	Offset	In/out mm	Width (M/D mm)	R/L	Non-convertible	
						0.022"	0.018"
Roth/MBT	-14°	+10°	0.5	4	UR	TDM2-62160e	TDM1-62160e
					UL	TDM2-62260e	TDM1-62260e
Roth	-25°	+4°	0.5	4	LL	TDR2-62360e	TDR1-62360e
					LR	TDR2-62460e	TDR1-62460e
MBT	-20°	0	0.5	4	LL	TDM2-62360e	TDM1-62360e
					LR	TDM2-62460e	TDM1-62460e
Edgewise	0	0	0.5	4	UR/LL	TDE2-62160e	TDE1-62160e
					UL/LR	TDE2-62260e	TDE1-62260e



BUCCAL TUBES

Product parameters (1st molar Double Combination Weldable) RS MEDICO · Product catalogue

1st Molar Double Combination Weldable									
Rx	Torque	Offset	In/out mm	Width (M/D mm)	R/L	Non-convertible		Convertible	
						0.022"	0.018"	0.022"	0.018"
Roth/MBT	-14°	+10°	0.5	4	UR	TDM2-63160	TDM1-63160	TDM2-53160	TDM1-53160
					UL	TDM2-63260	TDM1-63260	TDM2-53360	TDM1-53360
Roth	-25°	+4°	0.5	4	LL	TDR2-63360	TDR1-63360	TDR2-53360	TDR1-53360
					LR	TDR2-63460	TDR1-63460	TDR2-53460	TDR1-53460
MBT	-20°	0	0.5	4	LL	TDM2-63360	TDM1-63360	TDM2-53360	TDM1-53360
					LR	TDM2-63460	TDM1-63460	TDM2-53460	TDM1-53460
Edgewise	0	6°	0.5	4	UR/LL	TDE2-63160	TDE1-63160	TDE2-53160	TDE1-53160
					UL/LR	TDE2-63260	TDE1-63260	TDE2-53260	TDE1-53260



Product parameters (1st molar Triple weldable) RS MEDICO · Product catalogue

1st Molar Triple Weldable									
Rx	Torque	Offset	In/out mm	Width (M/D mm)	R/L	Non-convertible		Convertible	
						0.022"	0.018"	0.022"	0.018"
Roth/MBT	-14°	+10°	0.5	4	UR	TDM2-64160	TDM1-64160	TDM2-54160	TDM1-54160
					UL	TDM2-64260	TDM1-64260	TDM2-54360	TDM1-54360
Roth	-25°	+4°	0.5	4	LL	TDR2-64360	TDR1-64360	TDR2-54360	TDR1-54360
					LR	TDR2-64460	TDR1-64460	TDR2-54460	TDR1-54460
MBT	-20°	0	0.5	4	LL	TDM2-64360	TDM1-64360	TDM2-54360	TDM1-54360
					LR	TDM2-64460	TDM1-64460	TDM2-54460	TDM1-54460
Edgewise	0	0	0.5	4	UR/LL	TDE2-64160	TDE1-64160	TDE2-54160	TDE1-54160
					UL/LR	TDE2-64260	TDE1-64260	TDE2-54260	TDE1-54260



Buccal Tube (Bondable) in kit 4pcs/kit					
System	Kit Description	Convertible		Non-convertible	
		0.022"	0.018"	0.022"	0.018"
MBT	Upper Single Lower Single	TDM2-50160	TDM1-50160	TDM2-60160	TDM1-60160
	Upper Double Rectanglar Lower Single	TDM2-50260	TDM1-50260	TDM2-60260	TDM1-60260
	Upper One Round One Rectangular Lower Single	TDM2-50360	TDM1-50360	TDM2-60360	TDM1-60360
	Upper Triple Lower Single	TDM2-50460	TDM1-50460	TDM2-60460	TDM1-60460
	Upper Triple Lower Double Rectanglar	TDM2-50560	TDM1-50560	TDM2-60560	TDM1-60560
	Upper Triple Lower One Rectangular	TDM2-50660	TDM1-50660	TDM2-60660	TDM1-60660
	Upper Double Rectanglar Lower Double Rectanglar	TDM2-50760	TDM1-50760	TDM2-60760	TDM1-60760
	Upper Double Rectanglar Lower Double Rectanglar	TDM2-50760	TDM1-50760	TDM2-60760	TDM1-60760
Roth	Upper Single Lower Single	TDR2-50160	TDR1-50160	TDR2-60160	TDR1-60160
	Upper Double Rectanglar Lower Single	TDR2-50260	TDR1-50260	TDR2-60260	TDR1-60260
	Upper One Round One Rectangular Lower Single	TDR2-50360	TDR1-50360	TDR2-60360	TDR1-60360
	Upper Triple Lower Single	TDR2-50460	TDR1-50460	TDR2-60460	TDR1-60460
	Upper Triple Lower Double Rectanglar	TDR2-50560	TDR1-50560	TDR2-60560	TDR1-60560
	Upper Triple Lower One Rectangular	TDR2-50660	TDR1-50660	TDR2-60660	TDR1-60660
	Upper Double Rectanglar Lower Double Rectanglar	TDR2-50760	TDR1-50760	TDR2-60760	TDR1-60760
Edge-wise	Upper Single Lower Single	TDE2-50160	TDE1-50160	TDE2-60160	TDE1-60160
	Upper Double Rectanglar Lower Single	TDE2-50260	TDE1-50260	TDE2-60260	TDE1-60260
	Upper One Round One Rectangular Lower Single	TDE2-50360	TDE1-50360	TDE2-60360	TDE1-60360
	Upper Triple Lower Single	TDE2-50460	TDE1-50460	TDE2-60460	TDE1-60460
	Upper Triple Lower Double Rectanglar	TDE2-50560	TDE1-50560	TDE2-60560	TDE1-60560
	Upper Triple Lower One Rectangular	TDE2-50660	TDE1-50660	TDE2-60660	TDE1-60660
	Upper Double Rectanglar Lower Double Rectanglar	TDE2-50760	TDE1-50760	TDE2-60760	TDE1-60760

Product features

Molar Band



Bionic shape design

fit well

Sand-blasted inner surface

Better adhesion

Suitable hardness and softness

Good size consistency

Laser marked perimeter

Clear and permant



Product parameters (1st molar band "rough inner surface")

RS MEDICO · Product catalogue

1st molar band (rough inner surface)							
Perimeter mm	RS Size	GAC Size	16	26	36	46	Item in set/4pcs
31	31	3	D60021-31	D60022-31	D60023-31	D60024-31	D60020-31
31.5	31+	4	D60021-31+	D60022-31+	D60023-31+	D60024-31+	D60020-31+
32	32	5	D60021-32	D60022-32	D60023-32	D60024-32	D60020-32
32.5	32+	6	D60021-32+	D60022-32+	D60023-32+	D60024-32+	D60020-32+
33	33	7	D60021-33	D60022-33	D60023-33	D60024-33	D60020-33
33.5	33+	8	D60021-33+	D60022-33+	D60023-33+	D60024-33+	D60020-33+
34	34	9	D60021-34	D60022-34	D60023-34	D60024-34	D60020-34
34.5	34+	10	D60021-34+	D60022-34+	D60023-34+	D60024-34+	D60020-34+
35	35	11	D60021-35	D60022-35	D60023-35	D60024-35	D60020-35
35.5	35+	12	D60021-35+	D60022-35+	D60023-35+	D60024-35+	D60020-35+
36	36	13	D60021-36	D60022-36	D60023-36	D60024-36	D60020-36
36.5	36+	14	D60021-36+	D60022-36+	D60023-36+	D60024-36+	D60020-36+
37	37	15	D60021-37	D60022-37	D60023-37	D60024-37	D60020-37
37.5	37+	16	D60021-37+	D60022-37+	D60023-37+	D60024-37+	D60020-37+
38	38	17	D60021-38	D60022-38	D60023-38	D60024-38	D60020-38
38.5	38+	18	D60021-38+	D60022-38+	D60023-38+	D60024-38+	D60020-38+
39	39	19	D60021-39	D60022-39	D60023-39	D60024-39	D60020-39
39.5	39+	20	D60021-39+	D60022-39+	D60023-39+	D60024-39+	D60020-39+
40	40	21	D60021-40	D60022-40	D60023-40	D60024-40	D60020-40
40.5	40+	22	D60021-40+	D60022-40+	D60023-40+	D60024-40+	D60020-40+
41	41	23	D60021-41	D60022-41	D60023-41	D60024-41	D60020-41
41.5	41+	24	D60021-41+	D60022-41+	D60023-41+	D60024-41+	D60020-41+
42	42	25	D60021-42	D60022-42	D60023-42	D60024-42	D60020-42
42.5	42+	26	D60021-42+	D60022-42+	D60023-42+	D60024-42+	D60020-42+
43	43	27	D60021-43	D60022-43	D60023-43	D60024-43	D60020-43
43.5	43+	28	D60021-43+	D60022-43+	D60023-43+	D60024-43+	D60020-43+
44	44	29	D60021-44	D60022-44	D60023-44	D60024-44	D60020-44
44.5	44+	30	D60021-44+	D60022-44+	D60023-44+	D60024-44+	D60020-44+



Product parameters (2nd molar band"rough inner surface")

RS MEDICO · Product catalogue

2nd molar band (rough inner surface)							
Perimeter mm	RS Size	GAC Size	16	26	36	46	Item in set/4pcs
28	28	1	D70021-28	D70022-28	D70023-28	D70024-28	D70020-28
28.5	28+	2	D70021-28+	D70022-28+	D70023-28+	D70024-28+	D70020-28+
29	29	3	D70021-29	D70022-29	D70023-29	D70024-29	D70020-29
29.5	29+	4	D70021-29+	D70022-29+	D70023-29+	D70024-29+	D70020-29+
30	30	5	D70021-30	D70022-30	D70023-30	D70024-30	D70020-30
30.5	30+	6	D70021-30+	D70022-30+	D70023-30+	D70024-30+	D70020-30+
31	31	7	D70021-31	D70022-31	D70023-31	D70024-31	D70020-31
31.5	31+	8	D70021-31+	D70022-31+	D70023-31+	D70024-31+	D70020-31+
32	32	9	D70021-32	D70022-32	D70023-32	D70024-32	D70020-32
32.5	32+	10	D70021-32+	D70022-32+	D70023-32+	D70024-32+	D70020-32+
33	33	11	D70021-33	D70022-33	D70023-33	D70024-33	D70020-33
33.5	33+	12	D70021-33+	D70022-33+	D70023-33+	D70024-33+	D70020-33+
34	34	13	D70021-34	D70022-34	D70023-34	D70024-34	D70020-34
34.5	34+	14	D70021-34+	D70022-34+	D70023-34+	D70024-34+	D70020-34+
35	35	15	D70021-35	D70022-35	D70023-35	D70024-35	D70020-35
35.5	35+	16	D70021-35+	D70022-35+	D70023-35+	D70024-35+	D70020-35+
36	36	17	D70021-36	D70022-36	D70023-36	D70024-36	D70020-36
36.5	36+	18	D70021-36+	D70022-36+	D70023-36+	D70024-36+	D70020-36+
37	37	19	D70021-37	D70022-37	D70023-37	D70024-37	D70020-37
37.5	37+	20	D70021-37+	D70022-37+	D70023-37+	D70024-37+	D70020-37+
38	38	21	D70021-38	D70022-38	D70023-38	D70024-38	D70020-38
38.5	38+	22	D70021-38+	D70022-38+	D70023-38+	D70024-38+	D70020-38+
39	39	23	D70021-39	D70022-39	D70023-39	D70024-39	D70020-39
39.5	39+	24	D70021-39+	D70022-39+	D70023-39+	D70024-39+	D70020-39+
40	40	25	D70021-40	D70022-40	D70023-40	D70024-40	D70020-40
40.5	40+	26	D70021-40+	D70022-40+	D70023-40+	D70024-40+	D70020-40+
41	41	27	D70021-41	D70022-41	D70023-41	D70024-41	D70020-41

BUCCAL TUBES

Product parameters (Bicuspid band"rough inner surface")

RS MEDICO · Product catalogue

Bi-cuspid Bands																								
Size	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Ref.No.	
Upper (perimeter mm)	20	20.5	21	21.5	22	22.5	23	23.5	24	24.5	25	25.5	26	26.5	27	27.5	28	28.5	29	29.5	30	30.5	D50011-XX	
Lower (perimeter mm)	18.5	19	19.5	20	20.5	21	21.5	22	22.5	23	23.5	24	24.5	25	25.5	26	26.5	27	27.5	28	28.5	29	D50013-XX	



Product features

Molar Band Preweld with Buccal Tube



Product parameters

RS MEDICO · Product catalogue

1st molar band preweld with double buccal tube (one round one rectangular) in bulk									
System	Torque	Offset	In/out mm	Width (M/D mm)	Tooth	Convertible		Non-convertible	
						0.022"	0.018"	0.022"	0.018"
MBT	-14°	10°	0.5	4.0	16	D61021-**-0DM2d	D61021-**-0DM1d	D61023-**-0DM2n	D61023-**-0DM1n
	-14°	10°	0.5	4.0	26	D61022-**-0DM2d	D61022-**-0DM1d	D61024-**-0DM2n	D61024-**-0DM1n
Roth	-14°	10°	0.5	4.0	16	D61021-**-0DR2d	D61021-**-0DR1d	D61023-**-0DR2n	D61023-**-0DR1n
	-14°	10°	0.5	4.0	26	D61022-**-0DR2d	D61022-**-0DR1d	D61024-**-0DR2n	D61024-**-0DR1n
Edge-wise	0°	0°	0.5	4.0	16/26	D61021-**-0DE2d	D61021-**-0DE1d	D61021-**-0DE2n	D61021-**-0DE1n
	0°	0°	0.5	4.0	36/46	D61023-**-0DE2d	D61023-**-0DE1d	D61023-**-0DE2n	D61023-**-0DE1n



Product parameters

RS MEDICO · Product catalogue

1st molar band preweld with single buccal tube in bulk									
System	Torque	Offset	In/out mm	Width (M/D mm)	Tooth	Convertible		Non-convertible	
						0.022"	0.018"	0.022"	0.018"
MBT	-14°	10°	0.5	4.0	16	D61021-**-0DM2b	D61021-**-0DM1b	D61021-**-0DM2a	D61021-**-0DM1a
	-14°	10°	0.5	4.0	26	D61022-**-0DM2b	D61022-**-0DM1b	D61022-**-0DM2a	D61022-**-0DM1a
	-20°	0°	0.5	4.0	36	D61023-**-0DM2b	D61023-**-0DM1b	D61023-**-0DM2a	D61023-**-0DM1a
	-20°	0°	0.5	4.0	46	D61024-**-0DM2b	D61024-**-0DM1b	D61024-**-0DM2a	D61024-**-0DM1a
Roth	-14°	10°	0.5	4.0	16	D61021-**-0DR2b	D61021-**-0DR1b	D61021-**-0DR2a	D61021-**-0DR1a
	-14°	10°	0.5	4.0	26	D61022-**-0DR2b	D61022-**-0DR1b	D61022-**-0DR2a	D61022-**-0DR1a
	-25°	4°	0.5	4.0	36	D61023-**-0DR2b	D61023-**-0DR1b	D61023-**-0DR2a	D61023-**-0DR1a
	-25°	4°	0.5	4.0	46	D61024-**-0DR2b	D61024-**-0DR1b	D61024-**-0DR2a	D61024-**-0DR1a
Edge-wise	0°	0°	0.5	4.0	16/26	D61021-**-0DE2b	D61021-**-0DE1b	D61021-**-0DE2a	D61021-**-0DE1a
	0°	0°	0.5	4.0	36/46	D61023-**-0DE2b	D61023-**-0DE1b	D61023-**-0DE2a	D61023-**-0DE1a



Product parameters

RS MEDICO · Product catalogue

1st molar band preweld with double rectangular buccal tube in bulk									
System	Torque	Offset	In/out mm	Width (M/D mm)	Tooth	Convertible		Non-convertible	
						0.022"	0.018"	0.022"	0.018"
MBT	-14°	10°	0.5	4.0	16	D61021-**-0DM2c	D61021-**-0DM1c	D61021-**-0DM2o	D61021-**-0DM1o
	-14°	10°	0.5	4.0	26	D61022-**-0DM2c	D61022-**-0DM1c	D61022-**-0DM2o	D61022-**-0DM1o
	-20°	0°	0.5	4.0	36	D61023-**-0DM2c	D61023-**-0DM1c	D61023-**-0DM2o	D61023-**-0DM1o
	-20°	0°	0.5	4.0	46	D61024-**-0DM2c	D61024-**-0DM1c	D61024-**-0DM2o	D61024-**-0DM1o
Roth	-14°	10°	0.5	4.0	16	D61021-**-0DR2c	D61021-**-0DR1c	D61021-**-0DR2o	D61021-**-0DR1o
	-14°	10°	0.5	4.0	26	D61022-**-0DR2c	D61022-**-0DR1c	D61022-**-0DR2o	D61022-**-0DR1o
	-25°	4°	0.5	4.0	36	D61023-**-0DR2c	D61023-**-0DR1c	D61023-**-0DR2o	D61023-**-0DR1o
	-25°	4°	0.5	4.0	46	D61024-**-0DR2c	D61024-**-0DR1c	D61024-**-0DR2o	D61024-**-0DR1o
Edge-wise	0°	0°	0.5	4.0	16/26	D61021-**-0DE2c	D61021-**-0DE1c	D61021-**-0DE2o	D61021-**-0DE1o
	0°	0°	0.5	4.0	36/46	D61023-**-0DE2c	D61023-**-0DE1c	D61023-**-0DE2o	D61023-**-0DE1o



BUCCAL TUBES



Product parameters

RS MEDICO · Product catalogue

1st molar band preweld with triple buccal tube in bulk							
System	Torque	Offset	In/out mm	Width (M/D mm)	Tooth	Convertible	
						0.022"	0.018"
MBT	-14°	10°	0.5	4.0	16	D61021-**-0DM2e	D61021-**-0DM1e
	-14°	10°	0.5	4.0	26	D61022-**-0DM2e	D61022-**-0DM1e
Roth	-14°	10°	0.5	4.0	16	D61021-**-0DR2e	D61021-**-0DR1e
	-14°	10°	0.5	4.0	26	D61022-**-0DR2e	D61022-**-0DR1e
Edge-wise	0°	0°	0.5	4.0	26	D61021-**-0DE2e	D61021-**-0DE1e
	0°	0°	0.5	4.0	26	D61023-**-0DE2e	D61023-**-0DE1e



Lingual Attachment for Bands



Lingual Cleats

5.0 mm A212-45
6.5 mm A212-45L




Lingual Button with Round Base

A211-41



Lingual Button with Elliptical Base

A211-42



Lingual Sheath with hook

Monoblock

A219-45BL



Seating Lug

A21a-45




Lingual Sheath without hook

Monoblock

A219-45


Orthodontic Face Mask

Single Shaft




Orthodontic Face Mask

Double shaft



Face Mask



Face Mask with Silicon Pad & Head Cap (Single Shaft)

A51-22Lb

Face Mask with Silicon Pad (Single Shaft)

A51-20Lb

Face Mask with Silicon Pad & Head Cap (Double Shafts)

A52-22Ab

Face Mask with Silicon Pad (Double Shafts)

A52-20Ab



Face Mask with foam pad (single Shaft) blue large

A51-10Lb

Face Mask with foam pad (single Shaft) blue small

A51-10Sb

Face Mask with foam pad (single Shaft) red large

A51-10Lr

Face Mask with foam pad (single Shaft) red small

A51-10Sr

ACCESSORIES



Face Mask with foam pad (Double Shafts) blue large

A52-10Lb

Face Mask with foam pad (Double Shafts) blue small

A52-10Sb

Face Mask with foam pad (Double Shafts) red large

A52-10Lr

Face Mask with foam pad (Double Shafts) red small

A52-10Sr

High Pull System

High Pull System---headgear with curved J hook	Small Size	A53-31S
High Pull System---headgear with curved J hook	Medium Size	A53-31M
High Pull System---headgear with curved J hook	Large Size	A53-31L
High Pull System---headgear with curved J hook	Super Large	A53-31X



Mandibular pull system	Small Size	A53-41S
Mandibular pull system	Medium Size	A53-41M
Mandibular pull system	Large Size	A53-41L
Mandibular pull system	Super Large	A53-41X





Combination/Universal Pull System Small Size	A53-51S
Combination/Universal Pull System Medium Size	A53-51M
Combination/Universal Pull System Large Size	A53-51L
Combination/Universal Pull System Super Large	A53-51X



Safety Neck Pad without module short 140mm

A54-0S

Safety Neck Pad without module medium 170mm

A54-0M

Safety Neck Pad without module long 200mm

A54-0L

Safety Neck Pad with module short 140mm heavy force

A54-1S

Safety Neck Pad with module medium 170mm heavy force

A54-1M

Safety Neck Pad with module long 200mm heavy force

A54-1L

Safety Neck Pad with module short 140mm light force

A54-2S

Safety Neck Pad with module medium 170mm light force

A54-2M

Safety Neck Pad with module long 200mm light force

A54-2L

Face Bow




Face Bow Asymmetric
A32-01



Face Bow Symmetric
A31-01

Five Loop Torque Spring




A91-101	Torque Spring Maxillary small
A91-102	Torque Spring Maxillary medium
A91-103	Torque Spring Maxillary large
A91-201	Torque Spring mandibular small
A91-202	Torque Spring mandibular medium
A91-203	Torque Spring mandibular large

Stainless Steel Torque Spring



	Small	Medium	Large
0.016"x0.022"	A91-054	A91-055	A91-056
0.017"x0.025"	A91-064	A91-065	A91-066
0.018"x0.022"	A91-074	A91-075	A91-076
0.018"x0.025"	A91-084	A91-085	A91-086
0.019"x0.025"	A91-094	A91-095	A91-096


Molar distalizer



Niti Open Coil Spring 0.012" 180mm,with molar distalizer 8mm length Left side.	A92-123gL
Niti Open Coil Spring 0.012" 180mm,with molar distalizer 8mm length Right side.	A92-123gR
Niti Open Coil Spring 0.012" 180mm,with molar distalizer 10mm length Left side.	A92-123hL
Niti Open Coil Spring 0.012" 180mm,with molar distalizer 10mm length Right side	A92-123hR
Niti Open Coil Spring 0.012" 180mm,with molar distalizer 12mm length Left side	A92-123iL
Niti Open Coil Spring 0.012" 180mm,with molar distalizer 12mm length Right side	A92-123iR

Kit Cocle(2 pos/kit)	
Description	Number
Niti Open Coil Spring 0.012" 180mm,with 2 molar distalizer 8mm length Left/Right each.	A92-123gLR
Niti Open Coil Spring 0.012" 180mm,with 2 molar distalizer 10mm length Left/Right each.	A92-123hLR
Niti Open Coil Spring 0.012" 180mm,with 2 molar distalizer 12mm length Left/Right each.	A92-123iLR

Niti Open Coil Springs



Length 180mm
Wire Diameter


Open Spring 0.010" 180mm	A92-113
Open Spring 0.012" 180mm	A92-123

Niti Closed Coil Springs



Length	Wire Diameter	
	.010"	.012"
6mm	A94-110	A94-120
9mm	A94-111	A94-121
12mm	A94-112	A94-122

Niti Closed Coil Springs with eyelets



NiTi Closed Spring With Same Eyelet						
type	0.010"			0.012"		
	6mm	9mm	12mm	6mm	9mm	12mm
A-A	A94-210a	A94-211a	A94-212a	A94-220a	A94-221a	A94-222a

NiTi Closed Spring With Different Eyelet						
type	0.010"			0.012"		
	6mm	9mm	12mm	6mm	9mm	12mm
A-B	A94-310d	A94-311d	A94-312d	A94-320d	A94-321d	A94-322d
A-C	A94-310e	A94-311e	A94-312e	A94-320e	A94-321e	A94-322e

Crimpable Hook



Description	Length/mm	Width/mm	Slot 0.022"
Crimpable Hook (General)	1.3	2.2	A12-11
Mini Crimpable Hook (General)	1.3	1.9	A12-11M




Description	Length/mm	Width/mm	Slot 0.022"
Crimpable Hook (Medium)	2.2	2.2	A12-12



Description	Length/mm	Width/mm	Slot 0.022"
Crimpable Hook (Long)	3.0	2.2	A12-13



Description	Length/mm	Width/mm	Slot 0.022"
Crimpable Stops	/	2.2	A12-14
Crimpable Stops	/	1.9	A12-14M



Description

Length/mm

Width/mm

Slot 0.022"

Crimpable Hook Long Curved (Left)

6.9

2.5


A12-15L

Mini Crimpable Hook Curved (Left)

6.9

1.9

A12-15LM



Description

Length/mm

Width/mm

Slot 0.022"

Crimpable Hook Long Curved (Right)

6.9

2.5


A12-15R

Mini Crimpable Hook Curved (Right)

6.9

1.9

A12-15RM



Description

Length/mm

Width/mm


Slot 0.022"

Crimpable Hook Long Spiral (Left)

9.9

2.5

A12-16L



Description

Length/mm

Width/mm


Slot 0.022"

Crimpable Hook Long Spiral (Right)

9.9

2.5

A12-16R



Description

Length/mm

Width/mm

Slot 0.022"

Crimpable Hook BI-direction Shape

9.7

2.5

A12-17



Crimpable Mini Stops

0.015"

0.019"

2.0mm


A13-24M

A14-24M

3.0mm

A13-24


A14-24



T Type Pin
A10-51



Closed Crimpable Hook
A12-41



Crimpable Hook (T Hook)

0.022"

1.4mm

A12-01

1.6mm

A12-01 L



Crimpable Hook (Ball Hook)

0.022"

2.3mm

A12-01A

3.0mm

A12-01AL




Cross Tubes

0.022"

A12-34


0.018"

A11-34




Activity Archwire Stop

A12-24



Activity Stop with Straight Hook

A12-25



Activity Stop with Curved Hook

A12-25R

A12-25L

Lingual Button & Lingual Cleat



Description

Rectangle Base with Bi-direction
hook and slot .022"

Bondable

A226-33A



Description

Rectangle Base with Curved Hook (Right)

Bondable

A225-33Ra

Rectangle Base with Curved Hook (Left)

A225-33La



Description




Round Base with Curved Hook(Right)

Bondable

A225-31R

Round Base with Curved Hook(Left)

A225-31L

Base Shape	Bondable	Weldable	
Lingual Button with Round (Laser Sand-blasted)	A221-21	/	
Lingual Button with Elliptical (Laser Sand-blasted)	A221-22	/	
Lingual Button with Round	A221-31	A211-41	
Lingual Button with Elliptical	A221-32	A211-42	
Lingual Button with Triangle Base	A221-34BD	/	
Lingual Button with Round	A222-31	/	
Lingual Button with Elliptical	A222-32	/	
Lingual Button with Round Base with Flat Ring	A223-31	/	
Lingual Button with Elliptical Base with Flat Ring	A223-32	/	
Lingual Button with Round Base with Circle Ring	A224-31	/	



Lingual Button with hollow-out

A221-11

IA04-081

ACCESSORIES



Lingual Button Ceramic

Round Base

A231-51

Elliptical Base

A231-52




Ceramic Lingual Hook

A232-51




Ceramic Multi-Hook

A236-51



Ceramic Rectangular Tube
A239-51



Tongue Tammer
A227-33



Lingual Bite Plate (bondable)
A228-33



Lingual Button with Traction Chain

Base Shape

Round

Elliptical

S.S Color

A223-31C

A223-32C

Gold Coated

A223-31CG

A223-32CG

14K Gold

A223-31CK

A223-32CK



Withrertical eyelet

Base Shape

Round

Elliptical

S.S Color


A224-31C

A224-32C

Gold Coated

A224-31CG

A224-32CG




Lingual Sheath with hook
A219-45BL



Lingual Sheath without hook
A219-45

Lingual Retainer



Lingual Retainer Upper 38mm A61-138

Lingual Retainer Upper 40mm A61-140

Lingual Retainer Upper 42mm A61-142

Lingual Retainer Upper 45mm A61-145

Lingual Retainer Upper 47mm A61-147


Lingual Retainer Lower 29mm A61-229

Lingual Retainer Lower 31mm A61-231

Lingual Retainer Lower 33mm A61-233

Lingual Retainer Lower 35mm A61-235

Lingual Retainer Lower 37mm A61-237



Lingual Retainer Universal 1.2mm

A62-012

Lingual Retainer Universal 2.0mm

A62-020



Lingual Retainer 20mm

A63-320

Lingual Retainer 22mm

A63-322

Lingual Retainer 24mm

A63-324

Lingual Retainer 26mm

A63-326

Lingual Retainer 28mm

A63-328


Lingual Retainer 30mm

A63-330

Palatal Arch Bar	Distal Loop 37mm	A81-137
Palatal Arch Bar	Distal Loop 39mm	A81-139
Palatal Arch Bar	Distal Loop 41mm	A81-141
Palatal Arch Bar	Distal Loop 43mm	A81-143
Palatal Arch Bar	Distal Loop 45mm	A81-145
Palatal Arch Bar	Distal Loop 47mm	A81-147
Palatal Arch Bar	Distal Loop 49mm	A81-149
Palatal Arch Bar	Distal Loop 51mm	A81-151
Palatal Arch Bar	Distal Loop 53mm	A81-153
Palatal Arch Bar	Distal Loop 55mm	A81-155



Palatal Arch Bar	Mesial Loop 37mm	A81-237
Palatal Arch Bar	Mesial Loop 39mm	A81-239
Palatal Arch Bar	Mesial Loop 41mm	A81-241
Palatal Arch Bar	Mesial Loop 43mm	A81-243
Palatal Arch Bar	Mesial Loop 45mm	A81-245
Palatal Arch Bar	Mesial Loop 47mm	A81-247
Palatal Arch Bar	Mesial Loop 49mm	A81-249
Palatal Arch Bar	Mesial Loop 51mm	A81-251
Palatal Arch Bar	Mesial Loop 53mm	A81-253
Palatal Arch Bar	Mesial Loop 55mm	A81-255



Ligature Tie




A71-20M	Colorful Ligature Tie in Mickey
A71-20K	Colorful Ligature Tie in Kitty
A71-20F	Colorful Ligature Tie in Flower
A71-10	Ligature Tie Clear 40*26pcs
A71-20	Ligature Tie in Mixed Color 40*26pcs
A71-30	Key Shape Ligature Tie in Mixed Color

Color	Normal Shape	Key Shape
 #1	A71-201	A71-301
 #2	A71-202	A71-302
 #3	A71-203	A71-303
 #4	A71-204	A71-304
 #5	A71-205	A71-305
 #6	A71-206	A71-306
 #7	A71-207	A71-307
 #8	A71-208	A71-308
 #9	A71-209	A71-309
 #10	A71-210	A71-310
 #11	A71-211	A71-311
 #12	A71-212	A71-312
 #13	A71-213	A71-313
 #14	A71-214	A71-314
 #15	A71-215	A71-315
 #16	A71-216	A71-316
 #17	A71-217	A71-317
 #18	A71-218	A71-318
 #19	A71-219	A71-319
 #20	A71-220	A71-320
 #21	A71-221	A71-321
 #22	A71-222	A71-322
 #23	A71-223	A71-323
 #24	A71-224	A71-324

 #25	A71-225	A71-325
 #26	A71-226	A71-326
 #27	A71-227	A71-327
transparent		
 #28	A71-228	A71-328
 #29	A71-229	A71-329
 #30	A71-230	A71-330
 #31	A71-231	A71-331
 #32	A71-232	A71-332
 #33	A71-233	A71-333
 #34	A71-234	A71-334
 #35	A71-235	A71-335
 #36	A71-236	A71-336
 #37	A71-237	A71-337
 #38	A71-238	A71-338
 #39	A71-239	A71-339
 #40	A71-240	A71-340
 #41	A71-241	A71-341
 #42	A71-242	A71-342
 #43	A71-243	A71-343
 #44	A71-244	A71-344
 #45	A71-245	A71-345
 #46	A71-246	A71-346
 #47	A71-247	A71-347
 #48	A71-248	A71-348
 #49	A71-249	A71-349

Power Chain


Power Chain RS MEDICO · Product catalogue



Color	Long	Short	Closed
 #1	A72-201L	A72-201S	A72-201C
 #2	A72-202L	A72-202S	A72-202C
 #3	A72-203L	A72-203S	A72-203C
 #4	A72-204L	A72-204S	A72-204C
 #5	A72-205L	A72-205S	A72-205C
 #6	A72-206L	A72-206S	A72-206C
 #7	A72-207L	A72-207S	A72-207C
 #8	A72-208L	A72-208S	A72-208C
 #9	A72-209L	A72-209S	A72-209C
 #10	A72-210L	A72-210S	A72-210C
 #11	A72-211L	A72-211S	A72-211C
 #12	A72-212L	A72-212S	A72-212C
 #13	A72-213L	A72-213S	A72-213C
 #14	A72-214L	A72-214S	A72-214C
 #15	A72-215L	A72-215S	A72-215C
 #16	A72-216L	A72-216S	A72-216C
 #17	A72-217L	A72-217S	A72-217C
 #18	A72-218L	A72-218S	A72-218C
 #19	A72-219L	A72-219S	A72-219C
 #20	A72-220L	A72-220S	A72-220C
 #21	A72-221L	A72-221S	A72-221C
 #22	A72-222L	A72-222S	A72-222C
 #23	A72-223L	A72-223S	A72-223C
 #24	A72-224L	A72-224S	A72-224C


 #25	A72-225L	A72-225S	A72-225C
 #26	A72-226L	A72-226S	A72-226C
 #27	A72-227L	A72-227S	A72-227C
transparent			
 #28	A72-228L	A72-228S	A72-228C
 #29	A72-229L	A72-229S	A72-229C
 #30	A72-230L	A72-230S	A72-230C
 #31	A72-231L	A72-231S	A72-231C
 #32	A72-232L	A72-232S	A72-232C
 #33	A72-233L	A72-233S	A72-233C
 #34	A72-234L	A72-234S	A72-234C
 #35	A72-235L	A72-235S	A72-235C
 #36	A72-236L	A72-236S	A72-236C
 #37	A72-237L	A72-237S	A72-237C
 #38	A72-238L	A72-238S	A72-238C
 #39	A72-239L	A72-239S	A72-239C
 #40	A72-240L	A72-240S	A72-240C
 #41	A72-241L	A72-241S	A72-241C
 #42	A72-242L	A72-242S	A72-242C
 #43	A72-243L	A72-243S	A72-243C
 #44	A72-244L	A72-244S	A72-244C
 #45	A72-245L	A72-245S	A72-245C
 #46	A72-246L	A72-246S	A72-246C
 #47	A72-247L	A72-247S	A72-247C
 #48	A72-248L	A72-248S	A72-248C
 #49	A72-249L	A72-249S	A72-249C

O-Ring



Length/mm	1/8"- 3.2mm	3/16"- 4.8mm	1/4"- 6.4mm	5/16"- 8mm	3/8"- 9.5mm
Light 2.5OZ	A72-001	A72-002	A72-003	A72-004	A72-005
Medium 3.5OZ	A72-006	A72-007	A72-008	A72-009	A72-010
Heavy 4.5OZ	A72-011	A72-012	A72-013	A72-014	A72-015

Other Elastic



Separator


A74-20



Rotation Wedge





A75-10

1 Glass Mirror.



Dimension details (Mirror made from Optical Glass be made from)

RS MEDICO · Product catalog



I61-01 Glass Mirror Lingual (Zeiss Lenses Material)



I61-02 Glass Mirror Buccal (Zeiss Lenses Material)



I61-03 Glass Mirror Occlusal for Adult (Zeiss Lenses Material)



I61-04 Glass Mirror Super Large (Zeiss Lenses Material)

1.Glass	Optical grade glass imported from Japan
2.Coating	Chromium 70% reflectivity, six layers of uniform coating, constant reflection image to ensure the realization
3.Coating plane	Plating up to 99.95% and shooting angle up to 140 degrees
4.Membrane edge	Adopt \3000r angle cutter, and the dimensional error of high-speed CNC polishing film edge treatment is less than 0.1mm
5.sample	Scratch proof, fall resistant, can be immersed in high temperature and high pressure for 120 degrees


1 Glass Mirror

Dimension details (Glass Mirror coated Titanium)


RS MEDICO · Product catalog




I62-00 Glass Mirror in set 5pcs/set




I62-01 Glass Mirror Lingual




I62-02 Glass Mirror Buccal



I62-03 Glass Mirror Occlusal for Child



I62-04 Glass Mirror Occlusal for Adult



I62-05 Glass Mirror Super Large

Glass Mirror Holder with movable tip

RS MEDICO · Product catalog



I34-01B Mirror Holder in blue

I34-01R Mirror Holder in red

I34-01K Mirror Holder in gold

I34-01G Mirror Holder in grey

2 Bracket Positioning Gauge

**I12-00**

Double-ended Bracket Positioning Gauge
0.022" (2pc/pk)

**I12-01**

Double-ended Bracket Positioning Gauge
0.022" 2.0-2.5mm, 3.0-3.5mm

**I12-02**

Double-ended Bracket Positioning Gauge
0.022" 3.5-4.0mm; 4.5-5.0mm

**I13-01**

Crossed Bracket Positioning Gauge

3 Arch Turret

**I20-002**

Arch Turret 0° 5 Slot

4 Typodont




I40-01 Resin Teeth Model with S.S Joint



I42-01 Orthodontic Training Typodont Set



I41-01A Resin Teeth with steel root



I42-03 Typodont (Stainless Steel Frame)




I42-02 Wax model, Angle Class II/1

5 Force Gauge



I81-01 Force Gauge


6 Tweezers



I31-01 Bracket Tweezer Anterior



I31-02 Bracket Tweezer Posterior



I32-01 Buccal Tube Tweezer

7 Bracket Opener



I90-003 Self-ligating bracket open tool



I90-005 Self-ligating bracket open tool mini size

8 Retractor

I51-01 O type lip retractor




C type lip retractor

RS MEDICO · Product catalog

	white	blue	clear	
small	I52-11W	I52-11B	I52-11C	
medium	I52-12W	I52-12B	I52-12C	
large	I52-13W	I52-13B	I52-13C	

M Type Lip Retractor

RS MEDICO · Product catalog

	white	blue	
small	I57-11W	I57-11B	
medium	I57-12W	I57-12B	
large	I57-13W	I57-13B	

Long handle cheek retractor

RS MEDICO · Product catalog

small	medium	large	lip	
I55-11	I55-12	I55-13	I55-14	

Anterior Lip Retractor

RS MEDICO · Product catalog

I56-11	small	
I56-12	large	
I56-13	double end	

T type lip retractor

RS MEDICO · Product catalog

	blue	clear	
small	I54-11B	I54-11C	
large	I54-13B	I54-13C	

9 Orthodontic BlackBackground Photo Image Contrast



I71-11

Contrast with handle Anterior Small



I71-12

Contrast with handle Anterior Large



I71-13

Contrast with handle Occlusal Small

Orthodontic Plier packaging



Protection design

Protect orthodontic
Plier during
transportation



Transparent window

clearly seen without
open the package

➤ Ligature Cutter 0 degree P211-101



Hardness of tip:HRC 58-63
Diameter up to 0.30mm(0.012")

➤ Distal End Cutter P211-102



Hardness of tip:HRC 58-63
Arch wires up to 0.021"*0.025" (0.53mm*0.64mm)

➤ Distal End Cutter with flush end P211-103



Hardness of tip:HRC 58-63
Arch wires up to 0.021"*0.025" (0.53mm*0.64mm)

➤ Light Wire Pliers with Cutter P211-104



Make different fine loops
for arch wires and spring loops.
Tip with cutting edge for wires cutting.

➤ Heavy Cutter P211-105



Jaw hardness > HRC58-63.
The cutting edge is made of cemented carbide
Used for cutting Below 0.30 "/0. 76mm.
Various Australian wire, nickel titanium.
wire, stainless steel wire.

➤ Ligature Cutter 8 degree P211-106



Cut off the redundant arch wire, remove the
bonding bracket, band ring.

➤ Ligature Cutter 13 degree P211-107



It is used to cut soft tie wire, rubber band, etc
The blade is thin and sharp, and the clamp tip is finely ground.
Cut tie wire capacity: 0.012" / 0.30mm

➤ Filament forming pliers P211-108



Used for bending fine curves of various specifications of bow wires with different radii and bending various spring curves.
Bending arch wire can reach: 0.30" / 7.6mm.

➤ Torque Pliers P211-201



Used in pair to make torques on rectangle wires,
thickness of tips: 0.050" (1.27mm)

➤ Crimpable Hook Pliers P211-202



To fix the crimpable hooks along arch wires.

➤ Anterior Bracket Removing Pliers P211-203




Used to remove anterior brackets.

➤ Band Removing Pliers P211-204



Used to remove bands and other bondable attachment.

➤ Niti Wire Cinch Back Pliers P211-205



Used to cinch back distal end of arch wires to avoid any injury of soft tissue. Niti wires' diameter up to 0.021" (0.53mm)

➤ Three Jaw Pliers P211-206



To make V loops and bend arch wires. Wire's diameter up to 0.030" (0.76mm)

➤ Elastic Pliers P211-207



For expanding elastic separating rings before separating molar gap.

➤ Weingart Pliers **P211-208**



Used to cinch back distal end of arch wires to avoid any injury of soft tissue

➤ How's Pliers **P211-209**



To hold arch wires for guiding it into buccal tubes, to hold ligature wires, or to form anterior band

➤ Bondage Removing Pliers **P211-210**




Used to remove bondage on the crown without any injury on enamel.

➔ Tweed Pliers P211-211




To make small loops, distal stops and Ω loops on arch wires

➔ Kim Pliers P211-212



To bend multi-loop edgewise wires, with cutter on tips.

➔ Loop Bending Pliers P211-213



To make different small loops for multi-loop edgewise arch wires, diameter of loops : 2.00, 3.20, 4.20 mm

➤ Omega Pliers **P211-214**



To make fine Ω loops. Three different loops diameter 0.063", 0.075", 0.10" (diameter 1.6mm, 2.0mm, 2.6mm)

➤ Posterior Bracket Removing Pliers **P211-216**



To remove the posterior bracket.

➤ Buccal Tube Removing Pliers **P211-217**




To remove the caps of convertible tubes

➤ Nancy Loop Forming Pliers P211-218



To make vertical loops of different height along arch wires, 3.0, 4.0, 5.0, 6.0mm Wires' diameter up to: 0.022" (0.56mm)

➤ Tear Drop Pliers P211-219



Teardrop shaped notches can be clamped on the gingival edge of any transparent dental braces that require elastic traction. Teardrop shaped notches can hook the elastic skin ring well, making it easier for patients to use elastic traction. The angle and depth of teardrop shaped traction hook can be freely adjusted according to the direction of elastic traction. Also called traction hook hole shearing pliers.

➤ Horizontal Pressure Pliers P211-220




Realize the additional force of individual tooth root torque by clamping a dent on the appliance. Clip out the depression at the undercut or the gingival side of the attachment, and increase the retention of the transparent braces or retainers.

➔ Vertical Pressure Pliers **P211-221**



The vertical clamp can clamp the depression in the proximal or distal buccal/lingual side of the twisted tooth to solve the over correction of the twisted tooth. A pair of depressions can also be pressed at the relative positions of the buccal and lingual sides of the torsion teeth to form a torsion couple. You can also press a dent on a certain tooth position to generate a light thrust to correct slight lip and tongue misalignment.

➔ Hole Punch Pliers 6mm **P211-222**



A semicircular notch can be cut to provide space for bonding the lingual buckle on the tooth surface, so as not to affect the wearing of the dental braces.

It can cut off the parts that need soft tissue cushioning to prevent the gums from being compressed by the dental braces.

The semicircular notch can be placed on any tooth position where lingual buckle or bracket needs to be bonded.

➔ Hole Punch Pliers 3mm **P211-223**



A semicircular notch can be cut to provide space for bonding the lingual buckle on the tooth surface, so as not to affect the wearing of the dental braces.

It can cut off the parts that need soft tissue cushioning to prevent the gums from being compressed by the dental braces.

The semicircular notch can be placed on any tooth position where lingual buckle or bracket needs to be bonded.

➔ One Second Bending Pliers

P211-227 0.5mm

P211-228 0.75mm

P211-229 1.0mm





NITI Super-Elastic Arch Wires

NITI Super - Elastic Round Arch Wires (10pcs/pk)

RS MEDICO · Product catalogue

Round	Upper			Lower		
	Square	Oval	Taper	Square	Oval	Taper
.010"	W111-10U	W121-10U	W131-10U	W111-10L	W121-10L	W131-10L
.012"	W111-12U	W121-12U	W131-12U	W111-12L	W121-12L	W131-12L
.014"	W111-14U	W121-14U	W131-14U	W111-14L	W121-14L	W131-14L
.016"	W111-16U	W121-16U	W131-16U	W111-16L	W121-16L	W131-16L
.018"	W111-18U	W121-18U	W131-18U	W111-18L	W121-18L	W131-18L
.020"	W111-20U	W121-20U	W131-20U	W111-20L	W121-20L	W131-20L

NiTi Super-Elastic Rectangular Arch Wires (10pcs/pack)

RS MEDICO · Product catalogue

Rectangular	Upper			Lower		
	Square	Oval	Taper	Square	Oval	Taper
.016"x.016"	W111-1616U	W121-1616U	W131-1616U	W111-1616L	W121-1616L	W131-1616L
.016"x.022"	W111-1622U	W121-1622U	W131-1622U	W111-1622L	W121-1622L	W131-1622L
.017X.022"	W111-1722U	W121-1722U	W131-1722U	W111-1722L	W121-1722L	W131-1722L
.017X.025"	W111-1725U	W121-1725U	W131-1725U	W111-1725L	W121-1725L	W131-1725L
.018X.022"	W111-1822U	W121-1822U	W131-1822U	W111-1822L	W121-1822L	W131-1822L
.018X.025"	W111-1825U	W121-1825U	W131-1825U	W111-1825L	W121-1825L	W131-1825L
.019X.025"	W111-1925U	W121-1925U	W131-1925U	W111-1925L	W121-1925L	W131-1925L
.021X.025"	W111-2125U	W121-2125U	W131-2125U	W111-2125L	W121-2125L	W131-2125L

NITI Thermal Arch Wires

NiTi Thermal Round Arch Wires (5pcs/pack)

RS MEDICO · Product catalogue

Round	Upper			Lower		
	Square	Oval	Taper	Square	Oval	Taper
0.010"	W211-10U	W221-10U	W231-10U	W211-10L	W221-10L	W231-10L
0.012"	W211-12U	W221-12U	W231-12U	W211-12L	W221-12L	W231-12L
0.014"	W211-14U	W221-14U	W231-14U	W211-14L	W221-14L	W231-14L
0.016"	W211-16U	W221-16U	W231-16U	W211-16L	W221-16L	W231-16L
0.018"	W211-18U	W221-18U	W231-18U	W211-18L	W221-18L	W231-18L
0.020"	W211-20U	W221-20U	W231-20U	W211-20L	W221-20L	W231-20L

NiTi Thermal Rectangular Arch Wires (5pcs/pack)

RS MEDICO · Product catalogue

RECTANGULAR	Upper			Lower		
	Square	Oval	Taper	Square	Oval	Taper
.016"x.016"	W211-1616U	W221-1616U	W231-1616U	W211-1616L	W221-1616L	W231-1616L
.016"x.022"	W211-1622U	W221-1622U	W231-1622U	W211-1622L	W221-1622L	W231-1622L
.017X.022"	W211-1722U	W221-1722U	W231-1722U	W211-1722L	W221-1722L	W231-1722L
.017X.025"	W211-1725U	W221-1725U	W231-1725U	W211-1725L	W221-1725L	W231-1725L
.018X.022"	W211-1822U	W221-1822U	W231-1822U	W211-1822L	W221-1822L	W231-1822L
.018X.025"	W211-1825U	W221-1825U	W231-1825U	W211-1825L	W221-1825L	W231-1825L
.019X.025"	W211-1925U	W221-1925U	W231-1925U	W211-1925L	W221-1925L	W231-1925L
.021X.025"	W211-2125U	W221-2125U	W231-2125U	W211-2125L	W221-2125L	W231-2125L

NITI Reverse Arch Wires

NITI Reverse Round Arch Wires (2pcs/pk)

RS MEDICO · Product catalogue

ROUND	UPPER	LOWER
0.012"	W341-12U	W341-12L
0.014"	W341-14U	W341-14L
0.016"	W341-16U	W341-16L
0.018"	W341-18U	W341-18L
0.020"	W341-20U	W341-20L



NITI Reverse Rectangular Arch Wires (2pcs/pk)

RS MEDICO · Product catalogue

Rectangular	UPPER	LOWER
.016"x.016"	W341-1616U	W341-1616L
.016"x.022"	W341-1622U	W341-1622L
.017X.022"	W341-1722U	W341-1722L
.017X.025"	W341-1725U	W341-1725L
.018X.022"	W341-1822U	W341-1822L
.018X.025"	W341-1825U	W341-1825L
.019X.025"	W341-1925U	W341-1925L
.021X.025"	W341-2125U	W341-2125L



Stainless steel Arch Wires (10pcs/pk)

RS MEDICO · Product catalogue

Round	Oval		Natural	
	Upper	Lower	Upper	Lower
.012"	W402-12U	W402-12L	W403-12U	W403-12L
.014"	W402-14U	W402-14L	W403-14U	W403-14L
.016"	W402-16U	W402-16L	W403-16U	W403-16L
.018"	W402-18U	W402-18L	W403-18U	W403-18L
.020"	W402-20U	W402-20L	W403-20U	W403-20L
Rectangular				
.016"x .016"	W402-1616U	W402-1616L	W403-1616U	W403-1616L
.016"x .022"	W402-1622U	W402-1622L	W403-1622U	W403-1622L
.017"x .022"	W402-1722U	W402-1722L	W403-1722U	W403-1722L
.017"x .025"	W402-1725U	W402-1725L	W403-1725U	W403-1725L
.018"x .025"	W402-1825U	W402-1825L	W403-1825U	W403-1825L
.019"x .025"	W402-1925U	W402-1925L	W403-1925U	W403-1925L
.021"x .025"	W402-2125U	W402-2125L	W403-2125U	W403-2125L

Coated Supper-elastic Arch Wires (1pcs/pk)

RS MEDICO · Product catalogue

Round	Square		Oval		Taper	
	Upper	Lower	Upper	Lower	Upper	Lower
.012"	W501-12U	W501-12L	W502-12U	W502-12L	W503-12U	W503-12L
.014"	W501-14U	W501-14L	W502-14U	W502-14L	W503-14U	W503-14L
.016"	W501-16U	W501-16L	W502-16U	W502-16L	W503-16U	W503-16L
.018"	W501-18U	W501-18L	W502-18U	W502-18L	W503-18U	W503-18L
.020"	W501-20U	W501-20L	W502-20U	W502-20L	W503-20U	W503-20L

Cu-Alloy Niti Arch Wires (1pc/pk)

RS MEDICO · Product catalogue

Round	Oval		Natural		Damon*
	Upper	Lower	Upper	Lower	Upper/Lower
.013"	A-W-CU021301	A-W-CU021302	A-W-CU031301	A-W-CU031302	A-W-CU091301
.014"	A-W-CU021401	A-W-CU021402	A-W-CU031401	A-W-CU031402	A-W-CU091401
.016"	A-W-CU021601	A-W-CU021602	A-W-CU031601	A-W-CU031602	A-W-CU091601
.018"	A-W-CU021801	A-W-CU021802	A-W-CU031801	A-W-CU031802	A-W-CU091801
Rectangular	Upper	Lower	Natural	Natural	Upper/Lower
.014"x.025"	A-W-CU02142501	A-W-CU03142501	A-W-CU03142501	A-W-CU03142502	A-W-CU09142501
.016"x.025"	A-W-CU02162501	A-W-CU03162501	A-W-CU03162501	A-W-CU03162502	A-W-CU09162501
.018"x.025"	A-W-CU02182501	A-W-CU03182501	A-W-CU03182501	A-W-CU03182502	A-W-CU09182501

Beta Titanium Arch Wires (1pcs/pk)

RS MEDICO · Product catalogue

Round	Oval		Natural	
	Upper	Lower	Upper	Lower
.016"	A-W-ST021601	A-W-ST021602	A-W-ST031601	A-W-ST031602
.018"	A-W-ST021801	A-W-ST021802	A-W-ST031801	A-W-ST031802
.020"	A-W-ST022001	A-W-ST022002	A-W-ST032001	A-W-ST032002
Rectangular	Upper	Lower	Upper	Lower
.016"x.016"	A-W-ST02161601	A-W-ST02161602	A-W-ST03161601	A-W-ST03161602
.016"x.022"	A-W-ST02162201	A-W-ST02162202	A-W-ST03162201	A-W-ST03162202
.017"x.022"	A-W-ST02172201	A-W-ST02172202	A-W-ST03172201	A-W-ST03172202
.017"x.025"	A-W-ST02172501	A-W-ST02172502	A-W-ST03172501	A-W-ST03172502
.018"x.018"	A-W-ST02181801	A-W-ST02181802	A-W-ST03181801	A-W-ST03181802
.018"x.025"	A-W-ST02182501	A-W-ST02182502	A-W-ST03182501	A-W-ST03182502
.019"x.025"	A-W-ST02192501	A-W-ST02192502	A-W-ST03192501	A-W-ST03192502
.021"x.025"	A-W-ST02212501	A-W-ST02212502	A-W-ST03212501	A-W-ST03212502

Supper Elastic NiTi Dimpled Arch Wires




RS MEDICO · Product catalogue

Round	Oval		Natural		Damon*
	Upper	Lower	Upper	Lower	Upper/Lower
.012"	W110-12U	W110-12L	W120-12U	W120-12L	W150-12U
.014"	W110-14U	W110-14L	W120-14U	W120-14L	W150-14U
.016"	W110-16U	W110-16L	W120-16U	W120-16L	W150-16U
.018"	W110-18U	W110-18L	W120-18U	W120-18L	W150-18U
Rectangular	Upper	Lower	Natural	Natural	Upper/Lower
.014"x.025"	W110-1425U	W110-1425L	W120-1425U	W120-1425L	W150-1425U
.016"x.022"	W110-1622U	W110-1622L	W120-1622U	W120-1622L	W150-1622U
.016"x.025"	W110-1625U	W110-1625L	W120-1625U	W120-1625L	W150-1625U
.017"x.025"	W110-1725U	W110-1725L	W120-1725U	W120-1725L	W150-1725U
.018"x.025"	W110-1825U	W110-1825L	W120-1825U	W120-1825L	W150-1825U
.019"x.025"	W110-1925U	W110-1925L	W120-1925U	W120-1925L	W150-1925U




THERMAL NITI Dimpled Arch Wires

RS MEDICO · Product catalogue

Round	Square		Oval		Damon
	Upper	Lower	Upper	Lower	Upper/Lower
.013"	W210-12U	W210-12L	W220-12U	W220-12L	W250-12U
.014"	W210-14U	W210-14L	W220-14U	W220-14L	W250-14U
.016"	W210-16U	W210-16L	W220-16U	W220-16L	W250-16U
.018"	W210-18U	W210-18L	W220-18U	W220-18L	W250-18U
Rectangular	Upper	Lower	Natural	Natural	Upper/Lower
.014"x.025"	W210-1425U	W210-1425L	W220-1425U	W220-1425L	W250-1425U
.016"x.022"	W210-1622U	W210-1622L	W220-1622U	W220-1622L	W250-1622U
.016"x.025"	W210-1625U	W210-1625L	W220-1625U	W220-1625L	W250-1625U
.017"x.025"	W210-1725U	W210-1725L	W220-1725U	W220-1725L	W250-1725U
.018"x.025"	W210-1825U	W210-1825L	W220-1825U	W220-1825L	W250-1825U
.019"x.025"	W210-1925U	W210-1925L	W220-1925U	W220-1925L	W250-1925U



Preformed ligature wire	0.20mm	0.25mm	0.30mm	0.35mm
Long 155mm	A42-10L	A42-20L	A42-30L	A42-40L
Short 55mm	A42-10S	A42-20S	A42-30S	A42-40S
Short twisted 25mm	A42-10T	A42-20T	A42-30T	A42-40T



Kobayashi	0.30mm	0.35mm
Long 155mm	A43-30L	A43-40L
Short 55mm	A43-30S	A43-40S
Short twisted 25mm	A43-30T	A43-40T



RS MEDICO · Product catalogue

RS MEDICO

Professional Manufacturer Focusing on Ortho Bracket

Add:#6-1 Huapai Building,Liuxia,Xihu District,

Hangzhou,310023China

Web:www.rsmedico.com

<http://www.rsmedico.com>

2023.08 Version



CNC Milling Workmanship